



# clinell®

[www.clinell.com.au](http://www.clinell.com.au)

Near neutral  
pH with high  
material  
compatibility

## CLINELL UNIVERSAL RANGE

Low level instrument grade disinfectant wipes with dual disinfection and detergent action

Clinell Universal Wipes contain a patented formula with wide spectrum antimicrobial activity, effective against MRSA, VRE, *Acinetobacter*, *Klebsiella* (CRE), Norovirus, Influenza (H1N1) and Herpes Simplex Virus in one minute. A mix of biocides, with different mechanisms of action, designed to prevent bacterial resistance.

# UNIVERSAL WIPES

Clinell Universal Wipes provide a single universal product that cleans and disinfects in one easy step.



## RANGE OF USE

Suitable for use on surfaces and non-invasive medical devices. Broad compatibility with materials used throughout healthcare environments.

## TRUSTED BRAND

Developed by doctors, widely used in the UK and Australia. Over 164 million wipes have been sold in the last 12 months to Australian hospitals alone.

## ANTIMICROBIAL ACTIVITY


Low level instrument grade disinfectant wipes which act as both a disinfectant and detergent for surfaces. They contain a patented formula with wide spectrum antimicrobial activity, including bactericidal activity against MRSA, VRE, *Klebsiella* (CRE) and *Acinetobacter*, and virucidal activity against Norovirus, Influenza (H1N1) and Herpes Simplex Virus. A mix of biocides, with different mechanisms of action, designed to prevent resistance.

## SURFACE COMPATIBLE

Compatible with a wide range of surfaces and equipment. Approved by leading manufacturers for equipment and surface disinfection. Clinell Universal Wipes have a near neutral pH and are alcohol free. Visit: [clinell.com.au/compatibility](http://clinell.com.au/compatibility) for more information.

| BACTERIA                            | TIME    |
|-------------------------------------|---------|
| <i>Pseudomonas aeruginosa</i>       | 60 secs |
| <i>Enterococcus hirae</i>           | 60 secs |
| <i>Enterococcus faecium</i> (VRE)   | 60 secs |
| <i>Enterobacter cloacae</i> (CRE)   | 60 secs |
| <i>Klebsiella pneumoniae</i> (CRE)  | 60 secs |
| <i>Acinetobacter baumannii</i>      | 60 secs |
| <i>Staphylococcus aureus</i> (MRSA) | 60 secs |
| <i>Staphylococcus aureus</i>        | 60 secs |
| VIRUSES                             |         |
| Norovirus (MNV surrogate)           | 60 secs |
| Influenza (H1N1)                    | 60 secs |
| Herpes Simplex                      | 60 secs |

## How to clean surfaces effectively



Wear recommended PPE.



Remove one wipe from the pack. Apply to surface.



Working from top to bottom, clean to dirty, wipe in an 'S' shaped pattern, taking care not to go over the same area twice. Let the surface air dry for a minimum of 1 minute.




Change the wipe if it becomes dry or soiled and discard.

| PRODUCT                           | UNIT OF ISSUE | ORDER CODE   |
|-----------------------------------|---------------|--------------|
| Universal Wipes 40                | Pack of 40    | CW40AUS      |
| Universal Wipes 50 Clip Pack      | Pack of 50    | CWCP50AUS    |
| Universal Wipes 200               | Pack of 200   | CW200AUS     |
| Universal Wipes Tub 100           | Tub of 100    | CWTUB100AUS  |
| Universal Wipes Tub Refill 100    | Refill of 100 | CWTUB100RAUS |
| Universal Wipes Bucket 225        | Bucket of 225 | CWBUC225AUS  |
| Universal Wipes Bucket Refill 225 | Refill of 225 | CWBUC225RAUS |
| Wall Dispenser for CW200AUS       | Single unit   | CWD          |
| Wall Dispenser for CWTUB100AUS    | Single unit   | CWTUB100D    |
| Wall Dispenser for CWBUC225AUS    | Single unit   | CWBUC225D    |

USE DISINFECTANTS SAFELY.  
ALWAYS READ THE LABEL AND  
PRODUCT INFORMATION BEFORE USE.

ALWAYS FOLLOW MEDICAL DEVICE  
MANUFACTURER'S CLEANING AND  
DISINFECTING GUIDELINES.



ARTG No. 299351

GHA190002