

	Name and Style	S		M		L		XL		XL		Compression	Colour	
		Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe	Open Toe	Closed Toe			
						Full Calf	Full Calf			Full Calf	Full Calf			
DAILY	UltraSheer Stockings Knee High	75165-30	75165-25	75166-30	75166-25	75167-30	75167-25			75168-30	75168-25			● Natural Regular
			75408-02		75409-02		75410-02				75415-02			● Natural Petite
			75165-04		75166-04		75167-04				75168-04			● Black
			73633-00		73633-01		73633-02				73633-03			● Diamond Pattern Black
			75201-19		75202-19		75203-19				75204-19			● Natural Regular
			75435-02		75436-02		75437-02				75438-02			● Natural Petite
			75201-04		75202-04		75203-04				75204-04			● Black
			73634-00		73634-01		73634-02				73634-03			● Diamond Pattern Black
		75241-21		75242-21		75243-21				75244-21			● Natural	
		75241-04		75242-04		75243-04				75244-04			● Black	
		75169-22		75170-22		75171-22				75172-22			● Natural	
		75169-04		75170-04		75171-04				75172-04			● Black	
		73635-00		73635-01		73635-02				73635-03			● Diamond Pattern Black	
		75205-23		75206-23		75207-23				75208-23			● Natural	
		75205-04		75206-04		75207-04				75208-04			● Black	
		73636-00		73636-01		73636-02				73636-03			● Diamond Pattern Black	
		75245-20		75246-20		75247-20				75248-20			● Natural	
		75245-04		75246-04		75247-04				75248-04			● Black	
		77686-12		77686-13		77686-14				77686-15			● Natural	
		77686-04		77686-05		77686-06				77686-07			● Black	
		77687-12		77687-13		77687-14				77687-15			● Natural	
		77687-04		77687-05		77687-06				77687-07			● Black	
		77688-12		77688-13		77688-14				77688-15			● Natural	
		77688-04		77688-05		77688-06				77688-07			● Black	
	75173-12		75174-12		75175-12				75176-12			● Natural		
	75173-04		75174-04		75175-04				75176-04			● Black		
	73637-00		73637-01		73637-02				73637-03			● Diamond Pattern Black		
	75209-13		75210-13		75211-13				75212-13			● Natural		
	75209-04		75210-04		75211-04				75212-04			● Black		
	73638-00		73638-01		73638-02				73638-03			● Diamond Pattern Black		
	75249-13		75250-13		75251-13				75252-13			● Natural		
	75249-04		75250-04		75251-04				75252-04			● Black		
	75177-04		75178-04		75179-04				75180-04			● Natural		
	75177-01		75178-01		75179-01				75180-01			● Black		
	75213-06		75214-06		75215-06							● Natural		
	75213-03		75214-03		75215-03							● Black		
	75476-02		75477-02		75478-02				75479-02			● Black		
	75476-03		75477-03		75478-03				75479-03			● Khaki (Beige)		
	75484-02		75485-02		75486-02				75487-02			● Black		
	75484-03		75485-03		75486-03				75487-03			● Khaki (Beige)		
	75492-02		75493-02		75494-02				75495-02			● Black		
	75492-03		75493-03		75494-03				75495-03			● Khaki (Beige)		
	75253-01		75254-01		75255-01				75256-01			● Black		
	75257-01		75258-01		75259-01				75260-01			● Black		
	75265-01		75266-01		75267-01				75268-01			● Black		
	75285-00		75286-00		75287-00				75288-00			● Black		
	75261-00		75262-00		75263-00				75264-00			● Black		
	75269-00		75270-00		75271-00				75272-00			● Black		
	75147-03		75148-03		75149-03				75150-03			● Black		
	75147-02		75148-02		75149-02				75150-02			○ White		
	75155-03		75156-03		75157-03				75158-03			● Black		
	75155-02		75156-02		75157-02				75158-02			○ White		
	75411-03		75412-03		75413-03				75414-03			● Black		
	75411-02		75412-02		75413-02				75414-02			○ White		
	75289-80		75289-81		75289-82				75289-83			● Pink		
	75289-60		75289-61		75289-62				75289-63			● Royal blue		
	75339-10	75339-11	75340-10	75340-11	75341-10	75341-11	75343-10	75343-11	75342-10	75342-11	75344-10	75344-11	● Black	
	78050-04	78050-15	78051-04	78051-15	78052-04	78052-15	77975-15	77975-16	78053-04	78053-15	77976-15	77976-16	● Beige	
	78054-04	78054-15	78055-04	78055-15	78056-04	78056-15	77977-15	77977-16	78057-04	78057-15	77978-15	77978-16	● Beige	
	78042-08	78042-09	78043-08	78043-09	78044-08	78044-09			78045-08	78045-09			● Beige	
	78059-08	78059-09	78047-08	78047-09	78048-08	78048-09			78049-08	78049-09			● Beige	
	78476-10	78476-11	78477-10	78477-11	78478-10	78478-11			78479-10	78479-11			● Beige	
	78480-10		78481-10		78482-10				78483-10				● Beige	
	78031-05		78032-05		78033-05				78034-05				● Beige Left	
	78031-06		78032-06		78033-06				78034-06				● Beige Right	
	78035-05		78036-05		78037-05				78041-05				● Beige Left	
	78035-06		78036-06		78037-06				78041-06				● Beige Right	

Selecting the product that serves your patient's needs




Different types of graduated compression socks or stockings are required for different indications.

The **JOBST® Medical LegWear range** is designed to relieve tired, aching legs, help treat venous conditions, and help prevent swelling of the feet, ankles and legs.

A measured amount of graduated compression is applied to your legs. It is highest at the ankle and tapers off up the leg, promoting normal flow of venous blood up the leg. It is therefore ideal for patients who are mobile and requiring support of the venous system to reduce oedema.

Alternatively, anti-embolic stockings are intended for use during surgical procedures and on patients confined to bed. They are designed to prevent blood clots rather than active oedema management.



1 Select the required compression

15-20mmHg

- Tired aching legs
- Minor ankle, leg and foot swelling
- Minor varicose during pregnancy
- Long distance travel

20-30mmHg

- Tired aching legs
- Moderate to severe varicose during pregnancy
- Post vein surgery/post sclerotherapy
- Moderate tendency of oedema
- Superficial thrombophlebitis
- Helps prevent recurrence of ulcerations

30-40mmHg

- Severe varicose veins
- Lymphatic oedema
- Post vein surgery/post sclerotherapy
- Venous ulcers
- Manages manifestations of Post Thrombotic Syndrome (PTS), Chronic Venous Insufficiency CVI/CVD
- Aids in prevention of thrombotic syndrome

JOBST® Medical LegWear should not be worn and is contraindicated for the following conditions:

- Ischemia (eg. Advanced arterial disease)*
- Uncontrolled congestive heart failure
- Untreated septic phlebitis
- Phlegmasia coerulea dolens

*As part of a comprehensive assessment, a Doppler Ultrasound is used in the vascular screening process to clarify the vascular aetiology. A hand-held Doppler will assist in determining the degree, if any, of arterial disease. Using the Doppler the Ankle/Brachial Index (ABPI) can be calculated.

2 Select garment style



* open toe available


Some considerations when selecting the most appropriate style for your patient:

- Indication being treated/managed
- Lifestyle, mobility and age
- Ability to don and tolerate level of compression
- Skin Integrity
- Open toe or closed toe requirements
- Maternity Style only available in JOBST® UltraSheer

3 Measure for JOBST Medical LegWear

Measuring positions

- cH - Hip**
Circumference at largest dimension around buttocks.
- cG - Thigh**
Circumference at largest dimension.
- cC - Calf**
Circumference at largest dimension.
- cB - Ankle**
Circumference at smallest dimension.



Provide the following measurements for:

- Knee High - cB & cC
- Thigh High - cB, cC & cG
- Waist high - cB, cG & cH

- Each leg is measured independently
- Measure in the early morning when possible as swelling occurs later in the day

JOBST UltraSheer, forMen, forMen Casual, ActiveWear, Relief and Sport Sizing FC: Full Calf

Style	Measurement	Small	Medium	Large	Large FC	XLarge	XLarge FC
Unisex	Ankle	18 - 21cm	21 - 25cm	25 - 29cm	25 - 30cm	29 - 33cm	30 - 36cm
	Calf	28 - 38cm	30 - 42cm	32 - 46cm	46 - 61cm	34 - 50cm	46 - 61cm
	Thigh	40 - 62cm	46 - 70cm	54 - 78cm		60 - 81cm	
	Hip	71 - 117cm	76 - 127cm	81 - 137cm		102 - 166cm	
Maternity	Ankle	18 - 21cm	21 - 25cm	25 - 29cm		29 - 33cm	
	Thigh	40 - 62cm	46 - 70cm	54 - 78cm		60 - 81cm	
	Hip	81 - 132cm	86 - 142cm	91 - 152cm		112 - 180cm	

UltraSheer: Petite Knee High Length 34-39cm, Regular Knee High Length 39-45cm