

GAMMEX[®] Latex

S U R G I C A L G L O V E S

Latex, powder-free surgical glove offering unsurpassed combination of comfort, protection and sensitivity

KEY FEATURES AND BENEFITS:

- Proprietary glove formulation delivers a soft and sensitive yet robust surgical glove
- Made with Ansell's unique anatomic formers to enhance comfort and fit to help reduce hand fatigue
- Shaped finger pads improve fingertip sensitivity
- Additional cleaning processes ensure low residues of both latex proteins and chemical accelerators to minimize the risk of developing latex Type I and chemical Type IV allergies and skin sensitivities
- Tested for use with chemotherapy drugs and a wide range of chemicals used in hospitals*

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

**Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.*

Recommended For

- Damp donning
- Extra comfort during long, sustained procedures
- Tested for use with chemotherapy drugs
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES



GAMMEX[®] Innovation that Protects

Driven by innovation, GAMMEX[®] surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX[®] Latex

S U R G I C A L G L O V E S

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Straight with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Product Standards	AS/NZS 4179, ASTM D6978, ASTM F739, EN 16523-1, EN 374:2003 Microorganism, EN 374:2003 Low Chemical Resistant, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	<ul style="list-style-type: none"> • Xanthogen Polysulphides • Sodium di-n-butyldithiocarbamate (SDBC) • Zinc Diethyldithiocarbamate (ZDEC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9	9.5
Ansell Product Code	330048055	330048060	330048065	330048070	330048075	330048080	330048085	330048090	330048095

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 1-844-494-7854
 Email: insidesalesus@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu

Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21)2407 2286/88/62
 Email: infochina@ansell.com