



SOUTHERN CROSS



Supporting you

1300 362 923

www.healthtec.com.au

## Universal Examination Table

The SC electric lift Examination Table makes switching to an electric height adjustable table more affordable than ever before.

Offering a generous lifting capacity of 250kg, the motor is covered by an impressive 5 year warranty and moves from HI-LO within 16 seconds enabling the Practitioner to get the patient positioned quickly.

Each end section can be raised to provide a backrest for left and right patient examination and procedures.



### Features

Electric Height Adjustment with footswitch

X-Frame Lift action minimises pinch and crush points

Table can be moved similar to moving a wheelbarrow

Durable powder coated rust resistant frame

High quality anti-microbial vinyl & CFC free foam

5 year warranty on frame & electrics

Australian designed & manufactured

### Specifications

250kg Maximum Lifting Capacity (moving)

400kg Maximum Working Load (stationary)

540mm Minimum Height

860mm Maximum Height

710mm Table Width

1930mm Table Length

Manual backrest adjustment +70° (both end sections)

Footswitch Control



### The Healthtec Warranty

All Healthtec™ manufactured items are covered for defects in material or workmanship for 5 years from date of purchase. Original suppliers warranty applies.

TGA AR-00106892-6

© August 2013






Supporting you

www.healthtec.com.au

# PRICE SHEET

## Universal Examination Table

| Product Description  | CODE                   | Qty |
|--|------------------------|-----|
| <b>Product Description</b>   |                        |     |
| SC Universal Examination Table (710 wide)  | 56221                  |     |
| <b>Foam &amp; Vinyl Options</b>  |                        |     |
| 75mm Comfort Foam  | 6005                   |     |
| Non-Stock Vinyl Option   | 6018                   |     |
| *please refer to our colour chart for vinyl colour options   |                        |     |
| <b>Nose Hole Option</b>  |                        |     |
| Nose Hole in Backrest  | 6011-2-7               |     |
| Nose Hole Plug   | 6010                   |     |
| <b>Armrest Option</b>  |                        |     |
| Procedure Arm  | 5102                   |     |
| <b>Gynae Option</b>  |                        |     |
| Gynae Leg Supports   | 5114                   |     |
| <b>Table Protection</b>  |                        |     |
| Paper Towel Roll Holder  | 6017                   |     |
| Disposable Fitted Bed Sheet – 10 per pack  | ANWFIT                 |     |
| <b>Extended Warranty</b>   |                        |     |
| 10 years Frame & Electrics   | 7800                   |     |
| <b>Freight</b>   |                        |     |
| Freight charges are not included in our prices.<br>Please contact your Healthtec representative for freight costs and options. |                        |     |
| Shipping Dimensions  | 2050mm x 850mm x 650mm |     |
| Shipping Weight  | 80kg (approximately)   |     |



**FRECALL**  
**1300 362 923**  
 info@healthtec.com.au

Healthtec Pty Ltd ABN 94 224 518 379 TGA AR-00106892-6  
 © August 2013

Disclaimer: Due to our policy of constant product improvement, specifications may change without notice. Check with us for the latest specifications and pricing.