

Instant Heat Pack


Non-woven lining helps prevent the pack from sweating and can be placed directly on the skin. Does not require a compress holder or clumsy ties. Needs no refilling and simply fold to activate therapy.

Product	Catalogue No.	Packaging
Instant Hot Pack	601301	40/Box

- Advanced formulation keep un-activated before use
- Components and agents used are non-toxic
- Extra strong bonded seals to avoid leakages
- Easy to follow instructions printed on pack
- Rounded edges, no sharp tags
- Wrap in soft cloth for patient comfort or
- Apply directly to the skin under supervision.
- Simply fold to activate therapy
- TGA #288882

Medsupply Heat Pack – Superior Design and Quality Materials

	Medsupply Heat Pack	Standard Heat Pack
Pouch Thickness	0.11mm	0.08mm
White Inner Carton	Yes	No
Bottom Protection layer	Yes	No
Shipping Carton Thickness	≥750g/m ²	<700g/m ²
Shipping Carton Layers	5	<5
Materials Composition	Sodium acetate	Magnesium sulfate
Shelf Life	2 years	1 year
Tempreture Range	49±2 °C	55-65°C
Temperature Chart	More Stable	Less Stable



Medsupply Advantages

- Sodium Acetate is a more stable chemical
- Sodium Acetate has a longer shelf life
- Sodium Acetate prevents skin injury (<50 °C)
- Stronger pouches, less chance of bursting (0.11mm)
- Safe patient handling instruction on inner & outer cartons
- Stronger cartons ensure safely transported packs