

TREATMENT TABLES

Southern Cross Treatment Table

Three Section

The Australian-made, Southern Cross electric lift examination table makes switching to a height adjustable table more affordable than ever before.

A cost effective sturdy Australian made table built to last offering a generous lifting capacity of 250kg. The motor is covered by an impressive 5 year warranty, and moves quietly from low to high in just 18 seconds, enabling the practitioner to get the patient positioned quickly.

The adjustable head section with nose hole can be slightly lowered enabling the prone patient to fully relax the neck muscles giving optimum comfort during treatments.

The adjustable backrest enables the table to quickly be turned into a chair for when you need your clients in a seated position.



Specifications

Maximum Working Load:	250kg dynamic load (moving) 400kg static load (stationary)
Minimum Working Height:	540mm
Maximum Working Height:	860mm
Table Width:	610mm
Table Length:	1840mm
Hi-Lo Time:	18 seconds
Cushioning:	50mm Enduro Foam
Covering:	Vinyl Upholstery (with PreFixx protective finish)
Base/Castors:	X-Frame with Wheels
Footswitch:	Single Footswitch (Hi-Lo function)
Frame:	Powder Coated Rust Resistant, Steel Frame



The Healthtec Warranty


All Healthtec™ manufactured items are covered for defects in material or workmanship for 5 years from date of purchase. Original suppliers warranty applies.

PRICE SHEET

Southern Cross Treatment Table

Three Section

Description	Code	Qty
Products		
Southern Cross Treatment Table (610mm width)	56251	
Options and Accessories		
Foam and Vinyl Upgrade		
Non-Stock Vinyl Upgrade	6018	
75mm Comfort Foam Upgrade	6005	
Head Section Options		
Nose Hole and Plug	6010	
Arm Rest Options		
Procedure Arm	5102	
Table Protection		
Fitted Table Cover	7204	
Warranty Options		
Extended Warranty — 10 yrs motor, 10 yrs frame	7800	
Freight		
Freight is calculated by a Healthtec Representative when you place your order.		
Shipping Dimensions	2050mm x 750mm x 650mm	
Shipping Weight	80kg	



FREECALL
1300 362 923
 info@healthtec.com.au

Healthtec Pty Ltd ABN 94 224 518 379 TGA AR-00106892-6
 © September 2012

Disclaimer: Due to our policy of constant product improvement, specifications may change without notice. Check with us for the latest specifications and pricing.