

# FOCUS

on general healthcare



**Influenza & Travel Vaccines** P.4



**Melanoma March** P.6



**Cosmetics & Aesthetics** P.12



**Women's & Men's Health** P.16



**Diagnostic Equipment** P.18

## CONTACT YOUR PERSONAL ACCOUNT MANAGER



**Chris Mitchell**

Brisbane & Southside, Ipswich to  
Gatton & Mackay to Ayr  
Mobile: 0448 031 326

**Chris Wands**

Sunshine Coast, Wide Bay  
& Capricorn  
Mobile: 0408 061 221



**Jeff Osbourne**

Metro & Regional NSW  
Mobile: 0428 759 080

**Michelle Warnecke**

Brisbane North, Brisbane CBD &  
Far North Queensland  
Mobile: 0419 289 852

**Melissa Kneale**

Gold Coast, Toowoomba &  
Northern Rivers NSW  
Mobile: 0407 031 599



**Amaan Vaddadi**

Metro & Regional VIC  
Mobile: 0400 121 010



ITEM ATTRACTS  
EXTRA SHIPPING  
CHARGES

# BACK TO BASICS

## WOUND CARE WORKSHOP

**HERVEY BAY**  
FRIDAY 21 MARCH 2025



**SCAN to secure your ticket now!**  
or visit: [ssaustralia.com.au/events](http://ssaustralia.com.au/events)

See back page for more information!

## Quick Tips for Choosing a Wound Dressing

Dressing selection can be overwhelming and it's no wonder with thousands of different dressings available. With such a wide variety of products to choose from, it's important that you know the types of dressings available, when and when not to use them, and how to use them effectively.

Fortunately, like medications, dressings typically fall into generic categories or classes. Understanding the generic category of the product, their actions, indications, contraindications, advantages and disadvantages for use, will help with decision making.

When evaluating the performance of your dressings, ask yourself these questions:

- Does the dressing stay intact and remain in place for an extended wear time?
- Prevent leakage between dressing changes?
- Does it provide a moist wound healing environment?
- Protect the periwound skin?
- Reduce wound-associated pain?
- Reduce or control wound odour?
- Does the person with the wound find the dressing conformable and comfortable?
- Is the dressing easy to apply and remove without causing trauma to the wound or periwound skin?
- Are the dressings cost-effective (remembering that cheap isn't always cost-effective)?
- Supported by evidence and resources?



Dr Michelle Gibb  
M Nsg Sci (NP) M Wound Care BN CF



Continue reading this article on our website!  
Scan the code, or visit [ssaustralia.com.au/news-and-blog](http://ssaustralia.com.au/news-and-blog)



### Wounds in a Flash!

Enhance your knowledge and memory of wound types through these fun, new flash card games! Available for purchase via the Wound Specialist Services website



[FREE WEBINAR]

## LYMPHOEDEMA CARE IN VASCULAR DISEASE MANAGEMENT

TUESDAY 25 FEBRUARY 2025, 1PM - 2PM AEST



**REGISTER NOW**

Sponsored by: the science of healthy skin **essity**

[FREE WEBINAR]

## WOUND CARE DRESSINGS EXPLAINED

THURSDAY 27 MARCH 2025, 1PM - 2PM AEST



**COMING SOON**

Scan the QR code for more information about our events

Sponsored by: **convatec** **solventum**

[FREE WEBINAR]

## ACCREDITATION TIPS AND Q&A SESSION

TUESDAY 6 MAY 2025, 1PM - 2PM AEST



**COMING SOON**

Scan the QR code for more information about our events

# Influenza and Travel Vaccines

## 2025 Influenza - Price Breaks Available!

<p><b>Fluzone® High-Dose Quadrivalent</b> Pack 5</p> <p>1530615 <b>\$260</b> ex GST</p>	<p><b>FLUCELVAX QUAD</b> Pack 10</p> <p>1530186 <b>\$180</b> ex GST</p>
<p><b>FluQuadri®</b> Pack 5</p> <p>1530178 <b>\$50</b> ex GST</p>	<p><b>Afluria Quad</b> Pack 10</p> <p>1530151 <b>\$99</b> ex GST</p>

### Diphtheria & Tetanus

<p><b>ADT™ Booster</b> Pack 5</p> <p>1263409 <b>\$99<sup>75</sup></b></p>
---

### Diphtheria, Tetanus & Pertussis & Polio

<p><b>ADACEL® polio</b></p> <p>1269609 <b>\$64<sup>13</sup></b></p>
---

### Diphtheria, Tetanus & Pertussis

<p><b>Boostrix</b></p> <p>1190062 <b>\$33<sup>25</sup></b></p>
<p><b>boostrixripv</b></p> <p>1263436 <b>\$59<sup>38</sup></b></p>
<p><b>Adacel</b></p> <p>1458817 <b>\$31<sup>83</sup></b></p>

### Polio

<p><b>IPOL®</b></p> <p>1263338 <b>\$45<sup>13</sup></b></p>
---

### Measles & Mumps Rubella

<p><b>Priorix®</b></p> <p>1312251 <b>\$25<sup>65</sup></b></p>
--

### Typhoid

<p><b>TYPHIM VI™</b></p> <p>1265323 <b>\$64<sup>13</sup></b></p>
--

### Human Papillomavirus (HPV)

<p><b>GARDASIL.9</b></p> <p>1386905 <b>\$202<sup>35</sup></b></p>
---

### Shingles

<p><b>ZOSTAVAX®</b> Zoster Vaccine Live</p> <p>1302817 <b>\$189<sup>05</sup></b></p>
--

### Meningococcal A/C/W-135 & Y

<p><b>MenQuadfi™</b></p> <p>1475657 <b>\$64<sup>13</sup></b></p>
--

### Meningococcal ACWY Conjugate

<p><b>Nimenrix™</b></p> <p>1323415 <b>\$58<sup>90</sup></b></p>
---

### Meningococcal ACWY Conjugate

<p><b>BEXSERO</b></p> <p>1341472 <b>\$115<sup>81</sup></b></p>
--

### Japanese Encephalitis

<p><b>IMOJEV</b></p> <p>1305697 <b>\$269<sup>80</sup></b></p>
<p><b>JESPECT®</b></p> <p>1273791 <b>\$91<sup>63</sup></b></p>

### Yellow Fever

<p><b>STAMARIL®</b></p> <p>1270548 <b>\$106<sup>64</sup></b></p>
--

### Rabies

<p><b>VERORAB</b></p> <p>1317237 <b>\$124<sup>55</sup></b></p>
<p><b>Rabipur®</b></p> <p>1265270 <b>\$110<sup>20</sup></b></p>

### Varicella

<p><b>VARIVAX™</b></p> <p>1263300 <b>\$57<sup>76</sup></b></p>
--

<p><b>Varilrix</b></p> <p>1425218 <b>\$54<sup>15</sup></b></p>
--

### Cholera

<p><b>DUKORAL®</b> Single</p> <p>1263347 <b>\$56<sup>24</sup></b></p>
<p><b>DUKORAL®</b> Twin</p> <p>1269046 <b>\$112<sup>48</sup></b></p>

### Hepatitis A

<p><b>AVAXIM™</b></p> <p>1265314 <b>\$45<sup>60</sup></b></p>
<p><b>havrix®</b> Paediatric</p> <p>1150728 <b>\$33<sup>25</sup></b></p>
<p><b>havrix®</b> Adult</p> <p>1136827 <b>\$46<sup>08</sup></b></p>
<p><b>VAQTA®</b> Paediatric</p> <p>1155742 <b>\$33<sup>25</sup></b></p>
<p><b>VAQTA®</b> Adult</p> <p>1263374 <b>\$43<sup>65</sup></b></p>

### Hepatitis B

<p><b>engerixb</b> Paediatric</p> <p>1270127 <b>\$12<sup>35</sup></b></p>
<p><b>engerixb</b> Adult</p> <p>1049525 <b>\$24<sup>94</sup></b></p>
<p><b>H-B-VAX II</b></p> <p>1265708 <b>\$18<sup>05</sup></b></p>

### Hepatitis A & B

<p><b>twinrix®</b> Adult</p> <p>1152230 <b>\$61<sup>75</sup></b></p>
<p><b>twinrix®</b> Junior</p> <p>1158940 <b>\$46<sup>55</sup></b></p>

**Contact your Personal Account Manager Today!**



## WORLD IMMUNIZATION WEEK 24 - 30 April 2025

World Immunization Week, celebrated in the last week of April, aims to highlight the collective action needed and to promote the use of vaccines to protect people of all ages against disease.

## Vaccine Fridges

**NEW IN!**



### Nuline Pharmacy Refrigerator

This QCPP approved, refrigerator unit is ideal to be used in pharmacies and hospitals. The fan-forced cooling ensures a stable temperature throughout the whole unit protecting your inventory from temperature fluctuations and comes with a data logger included.



1532456	Nuline Pharmacy Refrigerator Glass Door	570L	Each	<b>\$3,249<sup>00</sup></b>
1533908	Nuline Pharmacy Refrigerator Glass Door	135L	Each	<b>\$1,799<sup>00</sup></b>
1533916	Nuline Pharmacy Refrigerator Glass Door	350L	Each	<b>\$2,799<sup>00</sup></b>



### Medisafe Plus Vaccine Fridges

- QCPP cold chain accredited and TGA approved
- Digital temperature control system
- Dual temperature probes
- AQ Box cloud temperature monitoring with SMS and email alerts
- Audible and visual alarms
- Adjustable shelves
- Lockable door
- Low energy consumption



1279211	Medisafe Plus 80 Fridge	3 Shelf	80L	Each	<b>\$2,699<sup>00</sup></b>
1446355	Medisafe Plus 650 Fridge	5 Shelf	610L	Each	<b>\$7,699<sup>00</sup></b>
1244495	Medisafe Plus 371 Fridge	6 Shelf	381L	Each	<b>\$3,999<sup>00</sup></b>
1532812	Medisafe Plus 420 Fridge	6 Shelf	422L	Each	<b>\$4,099<sup>00</sup></b>
1415808	Medisafe Plus 1000 Fridge	10 Shelf	1000L	Each	<b>\$8,875<sup>00</sup></b>



### Medisafe Fridge Data Loggers

1445928	Logtag Premium USB Data Logger w- Display	Each	<b>\$129<sup>00</sup></b>
1454420	Clever Starter Kit	Each	<b>\$289<sup>00</sup></b>
1385917	Medisafe Fridge AQ Box Cloud Based Temperature Monitoring Device	Each	<b>\$649<sup>00</sup></b>



### Medisafe Fridge UPS and Battery Packs

1436813	UPS and Battery Pack CSCT-2000	Each	<b>\$2,899<sup>00</sup></b>
1512609	Large UPS and Battery Pack CSCT-3000	Each	<b>\$3,099<sup>00</sup></b>
1487701	UPS Battery Pack	Each	<b>\$1,369<sup>00</sup></b>



# MELANOMA MARCH

Melanoma March celebrates Melanoma Institute Australia's major annual fundraising campaign. Each Melanoma March brings together melanoma patients, their families, and local communities to raise awareness and funds for research to reach our goal of zero deaths from melanoma.

For more information, visit their website: <https://melanomamarch.org.au/>

## Dermatoscopes

### Skin Closure



#### Daclon Nylon Sutures 3/8 RC

SMI is recognised as a world-wide supplier of the highest quality surgical sutures. The Daclon range of nylon sutures have high tensile strength and are great for use in general medical applications.

- Soft passage through the tissue without sawing, tissue drag and trauma
- Excellent and permanent tissue support
- Guaranteed seal and product sterility
- Good handling properties
- Excellent knot tie down
- Minimal tissue reaction
- Absence of capillarity



**Scan Me**  
To explore our full range of sutures

1488155	6/0 DS 12mm 75cm Black	B12	\$27 <sup>00</sup>
1488473	6/0 DS 16mm 75cm Black	B12	
1488449	5/0 DS 12mm 75cm Black	B12	
1488430	5/0 DS 16mm 75cm Black	B12	
1488414	5/0 DS 19mm 75cm Black	B12	
1488341	4/0 DS 19mm 75cm Blue	B12	
1488392	4/0 DS 24mm 75cm Blue	B12	
1488376	3/0 DS 24mm 75cm Blue	B12	
1488465	2/0 DS 26mm 75cm Blue	B12	
1488503	2/0 DS 26mm 45cm Blue	B12	

### Skin Preparation



#### Riodine Povidone Iodine Solution 10%

- Prevent infection in cuts and abrasions.
- Treatment of infected wounds and minor skin infections (boils, infected nails, chicken pox sores, tinea and other fungal and bacterial infections).

1312863	Riodine Povidone Iodine Solution 10%	15ml	Each	\$5 <sup>31</sup>
1228796	Riodine Povidone Iodine Solution 10%	100ml	Each	\$6 <sup>26</sup>
1228802	Riodine Povidone Iodine Solution 10%	500ml	Each	\$11 <sup>33</sup>



#### Heine Delta 30 PRO Bundle

The DELTA 30 PRO is the high-end dermatoscope from HEINE.

- Extremely vivid, spatial impressions of lesions
- Choose the light colour to suit your own preferences with Heine colour shift
- Sturdy and Durable
- Intuitive and ergonomic
- Digital documentation options

1512579	Heine Delta 30 PRO Bundle	Each	\$3,795 <sup>25</sup>
1512587	Heine Delta 30 PRO	Each	\$3,130 <sup>25</sup>
1445227	Heine Delta 30 Charging Station	Each	\$517 <sup>75</sup>
1512560	Heine Distance Working Ring	Each	\$217 <sup>55</sup>
1490257	Heine Universal Smartphone Adaptor DELTA 30	Each	\$170 <sup>05</sup>
1445200	Heine Small Contact Plate Delta 30	Each	\$403 <sup>75</sup>



#### Heine Deltaone Dermatoscope

The HEINE DELTAone Dermatoscope stands out as a remarkable combination of size, clarity, and advanced technology. This pocket-sized device offers an array of features that make it an indispensable tool for dermatologists and healthcare professionals.

- Compact design
- A true field of view
- Provides 10x magnification
- Tailored illumination
- Rechargeable battery

1445235	Heine Deltaone Dermatoscope Standard Edition Black	Each	\$1,854 <sup>16</sup>
1446266	Heine Deltaone White Edition	Each	\$1,960 <sup>14</sup>
1485156	Heine DELTAone RoseWhite	Each	\$1,524 <sup>94</sup>
1488805	Heine Universal Smartphone Adaptor DELTAone	Each	\$177 <sup>75</sup>
1446258	Heine Small Contact Plate for Deltaone	Each	\$446 <sup>64</sup>

### Lighting



#### W.A LED General Examination Light GS300

The GS 300 Exam Light provides a whiter, brighter light than halogen and consumes less energy for improved facility efficiencies.

1280562	With Table/Wall Mount	Each	\$675 <sup>10</sup>
1280554	With Mobile Stand	Each	\$749 <sup>00</sup>



#### W.A LED General Examination Light GS600

The Green Series 600 Minor Procedure Light features three, white LEDs in a wider head with integrated handle for easy maneuverability.

1280570	With Table/Wall Mount	Each	\$1,315 <sup>50</sup>
1280589	With Mobile Stand	Each	\$1,401 <sup>70</sup>



**Conmed Hyfrecator 2000 with Autoclavable Pencil**

The ConMed Hyfrecator 2000 is a low-power, high-frequency, high-voltage electrosurgical unit (ESU), a medical device that is used for various skin surgical procedures. It is often used for the removal of skin cancers, skin tags, red skin, sun spots, skin lesions, warts, moles, and other unwanted tissue. It allows healthcare professionals to perform procedures with minimal bleeding and scarring.

1373463 Each \$2,621<sup>86</sup>



1273657	Conmed Hyfrecator Broadbase Disposable Tips	Non-Sterile	Blunt	B100	\$104 <sup>86</sup>
1169971	Conmed Hyfrecator Disposable Sterile Tips	Sterile	Blunt	B50	\$106 <sup>59</sup>
1273648	Conmed Hyfrecator Pinpoint Disposable Tips	Non-Sterile	Sharp	B100	\$104 <sup>89</sup>
1182089	Conmed Hyfrecator Tips	Sterile	Sharp	B50	\$106 <sup>59</sup>



**Bovie DERM 942 HF Desiccator Complete Kit**

The DERM 942 and its patent pending hand piece are designed and engineered to be the most reliable and durable desiccator available today. The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat.

- Bipolar & monopolar modes
- Output power adjusts by 1/10 watt increments up to 10 watts, then by 1 watt increments up to 40 watts
- More accurate power with Universal Power Control 100-240 VAC
- Ground reference unit means no grounding pad is required
- Large illuminated digital displays
- Universal cable fittings to accept most accessory cable plugs
- Made in USA. Relied on by thousands of surgeons
- 4 Year warranty

1390155 Each \$2,123<sup>54</sup>



1259827	Bovie Dermal Tips Disposable	Non-Sterile	Blunt	B100	\$112 <sup>72</sup>
1295322	Bovie Dermal Tips Disposable	Sterile	Blunt	B50	\$194 <sup>70</sup>
1273209	Bovie Dermal Tips Disposable	Non-Sterile	Sharp	B100	\$112 <sup>72</sup>
1295314	Bovie Dermal Tips Disposable	Sterile	Sharp	B50	\$194 <sup>70</sup>



**Walmay Cryospray WM227 Liquid Nitrogen Unit w- 5 Tips & Case**

The Walmay Cryospray uses liquid nitrogen to treat various medical conditions, particularly in dermatology and some surgical procedures. The extreme cold causes localized damage or the destruction of tissue, the treated tissue will slough off or be absorbed by the body's natural healing processes.

- The flask is balanced for even pouring and due to a higher centre of gravity
- The unit is suitable for left and right-hand use due to a movable regulator
- 5 tips included - 1 curved tip and 4 various cone tips
- High-quality stainless steel and brass components are used
- The steel probe produces high tensile strength
- Comes in a convenient carry bag for safe transport
- Flask size is 350ml
- New ergonomically shaped collar design

1483668 Each \$695<sup>00</sup>



**Walmay Liquid Nitrogen Extractor**

Ensure safety and reliability when working with liquid nitrogen with the Walmay Liquid Nitrogen Extractor.

1484117 Each \$179<sup>00</sup>



**Walmay Liquid Nitrogen Dewar 20ltr**

Designed specifically for liquid nitrogen storage. Characterised by its small calibers and low static evaporation loss.

1484125 Each \$999<sup>00</sup>



**Dewar Trolley 10-50 Litres**

Stainless Steel Trolley for moving Dewars. Two levels for carrying Dewars ranging from 10 up to 50 Litres.

1505866 Each \$329<sup>00</sup>



**Take Steps to Save Lives from Melanoma**

Melanoma March is a family-friendly 4km walk which provides a wonderful opportunity for the community to come together to honour lives lost to melanoma and support those impacted by the disease.

If you want to be involved in the Melanoma March this year, you can find your local march and information about donating by scanning the QR Code.



# Surgical Instruments



## Forceps

1499874	Forceps Dressing Adson 12cm Sayco Sterile Disposable	Each	\$2 <sup>66</sup>
1450123	Forceps Tissue Adson Fine 1x2 Teeth 12cm Sayco Sterile Disposable	Each	\$2 <sup>80</sup>
1499793	Forceps Tissue Adson Micro 1x2 Teeth 12cm Sayco Sterile Disp	Each	\$3 <sup>14</sup>
1148630	Forceps Tissue Adson 1x2 Micro 12cm ARMO (Clinic)	Each	\$24 <sup>75</sup>
1499815	Forceps Artery Halstead Mosquito Curved 12.5cm Sayco St Disposable	Each	\$5 <sup>46</sup>




## Scissors

1499785	Dressing Sh/BI Straight 12.5cm Sayco Sterile Disposable	Each	\$4 <sup>37</sup>
1499769	Iris Straight 11.5cm Sayco Sterile Disposable	Each	\$3 <sup>56</sup>
1499777	Iris Curved 11.5cm Sayco Sterile Disposable	Each	\$3 <sup>56</sup>



## Scan Me

To explore our full range of surgical instruments



## Kai Skin Bopsy Punch

Sterile single-use, stainless steel instrument for accurate and gentle skin biopsies.

- Seamless stainless steel sharp cutting edge
- Ribbed handle for sure grip and control
- Sterile and ready to use
- Disposable
- 5-year shelf life

1034110	Kai Skin Biopsy Punch 2mm BP-20F	B20	\$71 <sup>00</sup>
1057590	Kai Skin Biopsy Punch 3mm	B20	\$71 <sup>00</sup>
1072022	Kai Skin Biopsy Punch 4mm	B20	\$71 <sup>00</sup>
1234336	Kai Skin Biopsy Punch 5mm	B20	\$71 <sup>00</sup>
1069949	Kai Skin Biopsy Punch 6mm	B20	\$71 <sup>00</sup>
1070378	Kai Skin Biopsy Punch 8mm BP-80F	B20	\$71 <sup>00</sup>

# Update on I.V. Fluids



The global shortage of I.V. fluids continues to be in effect. This situation continues to affect the supply chain for all Sodium Chloride and some Hartmans Solution lines from Baxter, Fresenius Kabi, and B Braun.

We are pleased to relay that there has been some temporary relief on supply issues for Hartmans Solution, and product lines from Baxter and B Braun are in stock and available as listed below.

The monthly allocation/review of stock remains in place for all Baxter Saline product lines. This measure is expected to continue until at least the end of June 2025 and potentially through to the end of the year. The temporary alternative supply line of saline fluids from Lavoisier is still available and continues to be allocated weekly.

When determining size and quantities to order, please keep in mind that currently the product with the most demand is 250ml Saline. Due to its widespread use for injectables, no stock has been available to wholesalers for some time. The second most saline product in demand is 500ml Saline.

All products list below are currently available; however, lead times and quantities may vary as the situation evolves.

If you have any questions or require assistance, please do not hesitate to reach out to our friendly customer support team on 1800 777 518.

## Saline

### Sodium Chloride



1529706	B Braun Sodium Chloride 0.9% I.V. 500ml	Each	\$5 <sup>15</sup>
1529714	B Braun Sodium Chloride 0.9% I.V. 1000ml <i>Purchase limits increased</i>	Each	\$6 <sup>10</sup>
1532251	Sodium Chloride 0.9% 100ml I.V. Bag Lavoisier	Each	\$9 <sup>45</sup>
1529242	Sodium Chloride 0.9% 100ml I.V. PVC Bag Bioluz	Each	\$9 <sup>45</sup>
1529234	Sodium Chloride 0.9% 250ml I.V. Bag Lavoisier	Each	\$12 <sup>50</sup>
1529226	Sodium Chloride 0.9% 500ml I.V. Bag Lavoisier	Each	\$15 <sup>75</sup>

## Hartmans Solution

### Hartmans Solution I.V.

The Hartmanns and modified HARTMANN's IV Infusion is indicated as a source of water and electrolytes. It is also used in patients as a source of bicarbonate in the treatment of mild to moderate metabolic acidosis associated with dehydration or associated with potassium deficiency.



1036044	Hartmans Solution 500ml I.V.	Each	\$2 <sup>68</sup>
1036035	Hartmans Solution 1000ml I.V.	Each	\$2 <sup>58</sup>
1529722	B.Braun Compounded Sodium Lactate EP I.V. 1000ml <i>Purchase limits increased</i>	Each	\$6 <sup>65</sup>



### Face Mask Earloop Lvl 2 Pink B50

- Shingle pleat design prevents fluid pooling and cross contamination
- Made from high-quality materials that provide reliable protection without compromising breathability
- Ultrasonically sealed, soft, stretchy earloops free of natural rubber latex prevent pulling and irritation
- The adjustable nosepiece ensures a comfortable, customised fit

1480979 B50 \$5<sup>50</sup>



### Micro-Touch NitraFree Accelerator Free Pink Nitrile Gloves

The MICRO-TOUCH® NITRAFREE™ Nitrile examination gloves are made from a durable and soft nitrile formulation manufactured without chemical accelerators, helping to protect users from Type IV chemical allergies.

- Tested for use with chemotherapy drugs
- Standard patient examinations
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients

1400207	Small	B100	\$10 <sup>21</sup>
1455613	Medium	B100	
1455621	Large	B100	
1455648	Extra Large	B100	



### Hammer Nitrile Powder Free Glove Black N/S

GloveOn Ridley continues to be the most popular choice of powder free latex glove in our range, offering unquantifiable aqueous latex proteins reading, and non-sensitising and non-irritant properties.

- Fingertip textured
- Fully textured
- Powder free
- Made with natural rubber latex
- Protein level = 50µg/g
- Ambidextrous
- Standard Cuff
- Natural colour

1338951	Small	B100	\$10 <sup>50</sup>
1338978	Medium	B100	
1338986	Large	B100	



### Disposable G-String Thin White

Suitable to wear for a long period of time during beauty treatments

- White with pink border
- Constructed from a soft, non-woven cotton fabric
- Elasticated waistband for superior comfort
- One size fits all
- Thin

1266358 P100 \$31<sup>30</sup>



### Terumo Unisharp Fine Gauge Short Needle 0.3ml 8mm

Fine gauge with a short needle for more comfortable objections

- Ultra-thin wall needle
- Ultra-sharp double bevel for easy skin penetration
- Easy to read markings on the clear barrel promoted accurate dosing
- Low dead space design minimises waste
- Latex and DEHP free
- Individually wrapped
- Sterile

1471171	31G	B100	\$25 <sup>37</sup>
1530410	32G	B100	\$29 <sup>48</sup>



### Imatec Nano Needle Stamp

The Auratio® device has the ability to deliver the desired amount of active ingredients just below the dermis at a constant depth, pain-free and over a large area. This is a major advantage over the conventional syringe. Medical grade stainless steel ensures the structural integrity of the micro-needle. These needles are gold-plated to ensure a high bio-compatibility with human skin, thus preventing irritations and allergic reactions.

- Auratio® is packaged in sterile individual blister packs.
- Available in 2 sizes: 0.5mm and 1.0mm.
- Single Use.
- Sterile with Ethylene Oxide (ETO).

1444778	Imatec Nano Needle Stamp 0.50mm	Each	\$51 <sup>16</sup>
1443259	Imatec Nano Needle Stamp 1.00mm	Each	\$44 <sup>51</sup>



### Couch Electric HiLo 3 Section Detach Arm Rests & Wings White

- Electric Hi Lo operation by foot switch
- Luxurious soft and comfortable head cushion
- Contoured comfortable table top
- Sleek design with rounded frame legs for safety
- Removeable padded wing armrests
- Padded prone arm rest
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

1461559 Each \$1,695<sup>00</sup>



SSS Australia are experts in the Cosmetics and Aesthetic product range. Reach out to your Personal Account Manager for tailored advice and recommendations on our full range of cosmetic products.

# Measurement



## Scales Bodisure LED Digital Display Personal Scale 180Kg

- Tempered safety glass
- LED display
- Weight measurement up to 180kg
- Unit conversion switch kg-lb

1397001 Each **\$35<sup>00</sup>**



## Seca Digital Floor Scale Large Display 150kg (not glass) 803

Seca Clara 803 electronic flat scales with high-quality two-component rubber surface.

- Capacity: 150kg
- Graduation: 100g
- Dimensions (mm): 316W x 37H x 326D

1278967 Each **\$149<sup>00</sup>**



## Scale Digital High Capacity 250kg M322600

An electronic floor scale with a capacity of 250kg with 50g graduations. It has BMI calculation, mother/child function, an auto off feature and tap to start. The perfect scale for daily use at hospitals, clinics & GP's.

1372157 Each **\$199<sup>00</sup>**



## Scale Seca Electronic Flat Floor 200kg x 100g (813)

Seca Robusta 813 electronic flat scales with very high capacity.

- Capacity: 200kg
- Graduation: 100g
- Dimensions (mm): 433W x 47H x 373D

1268915 Each **\$229<sup>00</sup>**



## Height Measure Retractable

- Retractable wall Height Measure
- 200 cm
- Black

1115228 Each **\$16<sup>34</sup>**



## Eye Charts

1010389	Eye Chart Universal Standard	6 Metre	Each	<b>\$24<sup>32</sup></b>
1010360	Eye Chart Direct	3 Metre	Each	<b>\$22<sup>80</sup></b>
1010351	Eye Chart Child Standard Animal	6 Metre	Each	<b>\$24<sup>32</sup></b>
1201952	Eye Chart Illiterate	3 Metre	Each	<b>\$22<sup>80</sup></b>

# MEDSTOCK<sup>®</sup> SILICONE RANGE PROTECTIVE ON SKIN EFFECTIVE ON WOUNDS

## Affordable Comfort with Medstock's Silicone Range

Designed to support healing in medical clinics, aged care, and at home. Medstock's silicone range offers cost-effective solutions for even the most fragile skin.



Effective exudate management for a variety of wounds



Soft, breathable and comfortable on the skin



Atraumatic on dressing removal



Silicone adhesive seals wound edges, reducing maceration risk



**Scan Me**  
For more information

# Tough on germs, gentle on surfaces

**Whiteley** <sup>TM</sup>  
Protecting lives and livelihoods since 1933



1035708	500ml Squeeze	Each	<b>\$9<sup>10</sup></b>
1035691	500ml Trigger Spray	Each	<b>\$8<sup>80</sup></b>
1035717	5 Litre	Each	<b>\$36<sup>76</sup></b>



Viraclean<sup>®</sup> is the result of years of intensive research into advanced cleaning and disinfecting technology. Viraclean<sup>®</sup> has been specially formulated to provide hospital grade, hard surface disinfectant power which is proven to kill a range of viruses including SARS-CoV-2 (COVID-19), bacteria and fungi.

Proudly Australian made you can rely on Viraclean to protect your patients, staff and your reputation.





## W.A Kleenspec Single Use LED Vaginal Specula

- Provides continuous output with a uniform spot
- Crisp white light
- Sealed LED light source and battery to reduce patient risk
- 100% Acrylic
- Smooth molded edges designed to provide maximum patient comfort



1331078	Extra-Small	B24	\$231 <sup>05</sup>
1331086	Small	B24	
1331108	Medium	B24	
1424408	Large	B14	



## IUD Insertion Kit Disposable Sterile

- Disposable kit containing:
- 1 x Cheron Tweezer – grips the sponge to clean the cervix
  - 1 x Pozzi Forceps – to hold the cervix
  - 1 x Scissor with littauer blades (length 22 cm) – to cut the IUD thread

1296922 B10 \$206<sup>20</sup>



## Vasectomy Forceps Pointed Curved 13cm Hifa (Theatre)

- Durable stainless steel construction
- 13 cm length for optimal control
- Lightweight

1293176 Each \$110<sup>72</sup>



## Forceps Vasectomy Curved 12cm ARMO (Clinic)

- Vasectomy Forceps are ratcheted forceps with sharp tips used to pierce the scrotal skin in vasectomy procedures.

1447513 Each \$49<sup>07</sup>



## Clamp Vasectomy 3mm Tip 14cm ARMO (Clinic)

- Vasectomy Towel Clamps are ratcheted forceps used to in non-suture vasectomy procedures.

1447521 Each \$58<sup>57</sup>



## Specimen Jar 70ml Yellow Top Labelled P5744SL

- Labelled to record patient details
- Durable container
- Yellow screw cap
- Leak-proof

1036474 Each \$0<sup>29</sup>



## IUD Hook 260mm Disposable

- The IUD hook is sterile and can be disposed easily.
- The length of the product is 260mm.

1295756 B20 \$189<sup>15</sup>



## IUD Flexible Conti Thread Retriever

- The PM Retrieval Hook for IUD is a sterile single use instrument. It can be easily manually curved according to the Torphology of the patient.

1430955 B10 \$117<sup>71</sup>



## Pregnancy TestPack Plus with OBC H.C.G. Urine Cassette Alere

The Alere™ TestPack +Plus hCG Combo Test is a rapid immunoassay for the qualitative detection of human Chorionic Gonadotrophin (hCG) in urine or serum for early detection of pregnancy.

1031155 B20 \$135<sup>00</sup>



## Pregnancy Test Surestep HCG Urine

The SureStep® hCG test is a simple and reliable rapid chromatographic immunoassay for the detection of human chorionic gonadotropin (hCG) in urine to aid in the early detection of pregnancy.

1416200 B25 \$19<sup>30</sup>



## Hysterometer Disposable (Uterine Sound)

- Individual Sterile pack.
- For dilation of the endovervix, and measurements of the uterine cavity.
- Visible graduations offer clear measurements.
- Soft distal part on a rigid handle.
- Available in 3 sizes

1295705	Blue	3.3mm	CH10	B25	\$120 <sup>08</sup>
1295691	Grey	4.3mm	CH12	B25	
1294725	Green	4.67mm	CH14	B25	



## Clinitek Status & Urine Chemistry Analyser

The CLINITEK Status®+ Analyzer features automatic checks (Auto-Checks\*). The analyzer automatically checks each test strip for humidity exposure, common sample interferences† and strip identification for Siemens test strips. Together, these provide improved clinical information.

- Enhances confidence in clinical decisions with (Auto-Checks) performed on every test
- Ensures consistent results interpretation with automatic timing and reading of every urine specimen
- Ready for connectivity to data management solutions (DMS), via a simple upgrade kit
- Offers a comprehensive point-of-care urinalysis menu for routine testing, detection of early kidney disease and hCG pregnancy testing
- Easy to use and convenient

1296175 Each \$916<sup>75</sup>



## MultistiX-10SG Strips 100

- Reliable testing for UTI, diabetes, and kidney disorders
- 10 tests for important markers on a single strip

1005899 Bt100 \$53<sup>29</sup>



## Thermal Paper Rolls for Clinitek Analyser

The CLINITEK Status®+ and DCA Vantage™ Analysers utilise the same printer paper.

1163181 B5 \$81<sup>99</sup>

# Diagnostic Equipment

## Spot Monitors



### W.A Spot Vital Signs 4400 with NIBP Suretemp Plug 6

- Records blood pressure in 15 seconds
- BP averaging helps ensure accurate hypertension detection
- 7" color touchscreen for a simplified workflow
- Options for documenting patient observations (like respiration rate and BMI)

1463756 Each **\$1,831<sup>00</sup>**



### W.A Spot Vital Signs 4400 Mobile Stand


1463640 Each **\$744<sup>00</sup>**



### W.A Spot Vital Signs 4400 Desk Stand

1463632 Each **\$341<sup>00</sup>**

## ECG



### W.A CP150 Interpretive 12 Lead Resting ECG Touch Screen

- Effortlessly input patient data
- Previews of ECG reports before printing
- Opt for adult and pediatric
- Portable battery operation
- Streamlined EMR connectivity
- The device holds up to 100 test results

1314270 Each **\$4,198<sup>00\*</sup>**



### W.A CP150 Hospital Cart W/Out Cable Arm/Shelf

1316516 Each **\$1,325<sup>00</sup>**



### W.A Mortara ELI 230 12-Lead Rest ECG Interp USB Host Flash

- Real-time preview of ECG data
- Post-acquisition review of acquired ECG
- Captures best 10 seconds of ECG data from memory
- Full size printouts
- Compact, portable form factor

1420208 Each **\$2,433<sup>00\*</sup>**



### AMEDTEC 12 Lead Resting ECG

- Complete integration with best practice, PDF importation with medical director
- Reliable pacemaker detection
- Banana plug leads for each electrode connection
- Faulty application is easily identified by graphic indicator
- Freeze and scroll back option during examination

1454110 USB Each **\$2,277<sup>00\*</sup>**

1454102 Bluetooth Each **\$2,669<sup>00\*</sup>**

## ABPM



### AMEDTEC ABPM ECO Pack Inc Lge/Med Cuff Belt/Holder USB Cable

- Stable 24h ABPM measurement
- Multi-day recording (up to 72 hours)
- Fast, convenient preparation of the recorder
- Patient-friendly monitoring

1454056 Each **\$2,828<sup>00</sup>**

# Diagnostic Equipment

## BPM



### W.A ProBP 2000 Digital BP Device and Flexiport Cuff Adult 11

- Blood pressure readings are captured on inflation in about 20 seconds
- Meets AAMI standards for blood pressure accuracy

1407821 Each **\$242<sup>00</sup>**



### Omron HEM7156T Plus Blood Pressure Monitor

- Cuff wrapping guide
- Body movement detection
- Irregular heartbeat detection
- Hypertension indicator
- Intellisense technology
- Bluetooth connectivity

1474758 Each **\$172<sup>00</sup>**

## ABPI Doppler



### Hadeco Smartdop XT6 Includes Temperature Probe and 8Mhz Vascular Probe

- Hadeco Smartdop will perform Ankle Brachial Index (ABI) testing and Toe Brachial Index (TBI) testing for greater accuracy in diabetic patients

1435868 Each **\$5,651<sup>00</sup>**



### Hadeco SmartDop XT 4 Cuff with Software/USB

- The Hadeco Smartdop XT is a fully automatic total vascular testing system for ABI and TBI studies. Simple one-button operation via PC for up to 6-site automatic blood pressure measurements.

1336037 Each **\$3,619<sup>00</sup>**

## Spirometers



### MIR Spirobank 2 Basic USB Spirometer

- User-friendly technology
- Icon-based navigation menu
- FVC and VC tests with main parameters always activated.
- Easy to use
- Low cost


1367463 Each **\$2,374<sup>00</sup>**



### MIR Spirobank 2 Smart Spirometer

- Multifunctional portable spirometer with oximetry option
- Stand alone mode
- PC via Bluetooth and USB
- Long lasting re-chargeable battery

1496204 Each **\$2,987<sup>00</sup>**



### Easy On Pc Spirometer

- Real time curves and Paediatric incentives
- Intuitive PC-based solution
- Automated user guidance
- Z-score, LLN and % predicted for fast interpretation of results

1418505 Each **\$2,992<sup>00</sup>**



### Scout Spirometry Tube Wrapped Standard

1205530 B50 **\$156<sup>00\*</sup>**



### Vitalograph Bacterial Viral Filter For Spirolab/bank Latex F

1252370 B50 **\$188<sup>00</sup>**



sssaustralia HEALTHCARE SUPPLIES



WOUND SPECIALIST Services

# BACK TO BASICS

## WOUND CARE WORKSHOP

**HERVEY BAY**  
FRIDAY 21 MARCH 2025



**SCAN to secure your ticket now!**  
or visit: [sssaustralia.com.au/events](https://sssaustralia.com.au/events)

### IPSWICH

FRIDAY 2 MAY 2025

### TOWNSVILLE

FRIDAY 16 MAY 2025

### BALLINA

FRIDAY 6 JUNE 2025

### MACKAY

FRIDAY 27 JUNE 2025

### NORTHERN SYDNEY CBD

FRIDAY 25 JULY 2025

### ALBURY WODONGA

FRIDAY 22 AUGUST 2025

### BRISBANE CBD

FRIDAY 12 SEPTEMBER 2025

### ABOUT THE WORKSHOP

Focusing on the basics of wound care, Dr. Gibb will be discussing the principles of dressing selection, the challenges faced when treating wounds and reviewing case studies.

There will be several practical workshop activities and attendees will be able to use the time spent in the seminar towards their Continuing Professional Development.

### LEARNING OUTCOMES

Challenges in wound care  
Principles of dressing selection  
Workshop activities:

- Wound measurement and photography
- Barrier preparations and emollients
- Wound cleansing and biofilm
- Pressure injuries - staging
- Pressure off-loading & redistribution

Case study activity & discussion



### EDUCATOR

**Dr Michelle Gibb**

*M Nsg Sci (NP) M Wound Care BN CF*  
Director and Founder of  
Wound Specialist Services

**\$25 TICKETS**

**Includes**

- Lunch, morning tea, tea & coffee
- Digital certificate of attendance
- Goodie bag



sssaustralia HEALTHCARE SUPPLIES

1800 777 518

sales@sssaustralia.com.au

sssaustralia.com.au

Pricing Valid February 1, 2025 - May 31, 2025. All orders may be subject to a shipping fee. Pictures for illustration purposes only. All items while stocks last. All prices shown are exclusive of GST, unless marked\*. E&O.E.