



**SSS AUSTRALIA**  
HEALTHCARE SUPPLIES

## Aged Care 2022

# FURNITURE & EQUIPMENT CATALOGUE

SSS Australia aims to be your partner in your business. Whether it be a new site in your portfolio, your first venture, or upgrades to your existing furniture and equipment, we are here to help.

With experienced people in all areas of the process, from product selection to logistics, you can be assured we'll have the answers to your questions, throughout what can be a daunting process.

Talk to us today about your goals, and we'll be your partners in achieving them.

 **1800 777 518**

 **sales@sssaustralia.com.au**

 **www.sssaustralia.com.au**

# **SSS AUSTRALIA**

## **Supporting Quality Care**

SSS Australia is committed to supporting you in providing the highest quality of care for your residents. Our extensive range of furniture and equipment ensures we have the product solution for your facility. We choose to supply high quality, innovative products to help you provide the highest quality care for your residents.

For more information on any of the products in this booklet please do not hesitate to contact your personal account manager or our expert customer support team via phone 1800 777 518 or email- [sales@sssaustralia.com.au](mailto:sales@sssaustralia.com.au)

To assist in complying with your facility's COVID safe measures, we can also organise a digital meeting with you via Microsoft Teams. Talk to your personal account manager today.

### **CONTACT YOUR PERSONAL ACCOUNT MANAGER TO REQUEST A QUOTE:**

**Carol Edmond**  
**0408 181 938**  
**[carol.edmond@sssaustralia.com.au](mailto:carol.edmond@sssaustralia.com.au)**

**Fiona Carey**  
**0427 361 761**  
**[fiona.carey@sssaustralia.com.au](mailto:fiona.carey@sssaustralia.com.au)**

**Pam Harris**  
**0419 172 150**  
**[pam.harris@sssaustralia.com.au](mailto:pam.harris@sssaustralia.com.au)**



# TABLE OF CONTENTS

**DAILY LIVING AIDS**

**MANUAL HANDLING & SAFETY**

**MEDICATION & VACCINE FRIDGES**

**ARM CHAIRS**

**PRESSURE MANAGEMENT**

**LINEN & LAUNDRY**

**MEDICATION ADMINISTRATION**

**DIAGNOSTIC MEASURE**

**RESUSCITATION & CRITICAL CARE**

**PRICE  
MATCH  
POLICY**

SSS Australia is committed to remaining competitive on a comprehensive range of quality products from trusted brands. Our **Purchase With Confidence** policy means that we will make every effort to match genuine advertised prices from local competitors on the same stocked item.



## BATHROOM AIDS

### Aspire Shower Commodes

Easy-slide weight bearing platform footrest slides completely away for transfers.

- 8° backrest angle ensures patient comfort
- Improved toilet clearance
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- 4x total locking stainless steel castors provide brake access from all angles
- Safety arms prevent falls and swing back out of the way for transfers
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility




Shower Commode - Weight Bearing Platform, Att Propelled				
Code	Max Weight	Seat Height	Seat Width	Overall Width
1376160	200kg	590mm	460mm	535mm
1375989	200kg	590mm	530mm	605mm
1335162	400kg	590mm	600mm	725mm

Aspire Padded Commode Seats		
Code	Description	Fits
1376268	Grey Open Front 460mm	1376160, 1368540
1375997	Grey Open Front 530mm	1375989
1335170	Grey Open Front 600mm	1335162
1396250	Bowl and Pan Carrier	1376160,1368540, 1375989,1335162

### Padded Commode Seats

- Moulded and seamless soft PU seats designed for comfort and durability
- Easy clean, fully sealed surface for infection control
- Choose the size appropriately to fit your Aspire Shower Commode
- Adaptable to most shower commode brands



### Smik Shower Commode Sliding Footrest & Open Seat

- Can be positioned over the toilet for toileting
- Swingback arms assist side transfers
- Pan carrier included
- Padded seamless commode seat included
- Patient locking arms
- Weight bearing sliding footrest
- Slide in foot plate
- PVC slatted backrest

Code	Description
1454854	Smik Commode 550W
1386867	Smik Commode 450W



### Aspire Tilt-in-Space Shower Commode

- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Height adjustable, swing away footrests with calf pads for optimal positioning, support and ease of transfers
- Tool-free height adjustment to 4 settings allow the commode to be easily adapted to varying needs



Shower Commode - Weight Bearing Platform, Att Propelled				
Code	Max Weight	Seat Height	Seat Width	Overall Width
1368540	250kg	460-630mm	460mm	560mm

Shower Commode -Accessories	
1376268	Aspire Shower Commode Seat Open Front 460mm
1404857	Aspire Shower Commode Seat Closed Front 460mm

- Requires Seat Refer To Aspire Padded Commode



## Alphacare Lopital Reflex Electric Shower Commode

- Electrically operated height-adjustable settings (variable)
- Electrically operated tilt mechanism
- Easy-to-use touch controls on backrest
- Height-adjustable headrest
- Ergonomically designed backrest and seat.
- Height-adjustable legrests and calf supports
- Foldable swivel armrests with front over safety locking arms
- Suitable for use in the shower and toilet
- Carrier for stainless steel commode pan
- 4x braking castors
- Easy to recharge with battery charger provided
- Also available with Hand-control



Alphacare Lopital Reflex Shower Commode				
Code	Max Weight	Seat Height	Seat Width	Overall Width
1435310	135kg	490-1001mm	480mm	660mm

## Aquarius Mobile Shower Trolley

- Economical hydraulic trolley with 180kg SWL
- Powder coated stainless steel dropsides and endrails
- 795 W x 1950mm Long overall size including bumpers
- Liner inner (min) base size - 560 W x 1720 L x 180mm Deep
- Adjustable end rails with 120mm extra length
- Standard hose length; 1100mm, custom lengths available
- Height Range - 680 to 880mm
- Quality 125mm twin wheel castors with total lock and steer
- 135mm lowest clearance under the base
- Made in Australia



Code	Max Weight	Overall Size	Liner Base	Hose Length
1373250	180kg	795x1950mm	1720 x 560mm	1100mm

## Shower Chair Plastic Adjust Height Standard 120kg Aluminium

The Homecraft Aluminium Shower Chair enables comfortable & safe showering, with a curved seat for comfort & drainage. Made with a rust resistant aluminium frame. The uncoated internal sliding legs provide smooth effortless height adjustment, while avoiding any chance of water filling the legs. Reinforced frame using 25.4mm tubing enables users to gain a firmer grip making transferring safer and easier.

Code	Max Weight	Seat Height	Seat Width	Overall Width
1292633	110kg	390 - 540mm	440mm	590mm



## Shower Chairs

- Comfortable and hygienic plastic moulded seat
- Gentle slope allows water to easily flow away



Shower Chairs				
Code	Max Weight	Seat Height	Seat Width	Overall Width
1064094	125kg	420-550mm	440mm	570mm
1292641	125kg	420-570mm	440mm	645mm
1426508	200kg	375-550mm	440mm	512mm



## Over the Toilet Shower Chairs

Over the Toilet Shower Chairs				
Code	Max Weight	Seat Height	Seat Width	Overall Width
1175280	125kg	440-490mm	370mm	450mm
1214807	135kg	470-590mm	500mm	600mm
1214790	210kg	515mm	650mm	610mm



## Bariatric All-In-One Over Toilet & Shower Chair Commode

- Comes with square bowl with easy front access to remove bowl
- The seat lifts up and comes with a seat cover.
- Removable backrest adds stability when used as a bedside commode

Over the Toilet Shower Chairs				
Code	Max Weight	Seat Height	Seat Depth	Overall Width
1286250	295kg	465-570mm	450mm	790mm



## Aspire Deluxe Bedside Commode

- Allows for toileting in the bed room for users unable to travel between the bedroom and bathroom
- Available in Champagne or Slate Vinyl
- Easy Clean Vinyl Surfaces
- Removable Bowl with Lid & Carry Handle
- Max Weight: 160kg, Seat Height : 410-565mm,
- Seat Depth: 450mm, Seat Width: 460mm

Deluxe Bedside Commode	
Code	Colour
1408364	Slate (Blue)
1434810	Champagne



## Toilet Seat Reducer

- Designed for thin frail and bony patients who may slip down into the toilet bowl.
- Use the Seat Reducer on commode chairs, shower chairs or straight on a toilet seat.
- Soft and comfortable foam that will reduce the gap of the hole in the toilet seat.

Seat Reducer	
Code	Hole Size
1289616	20.5x8cm

## Pelican Shower Chair Cape

- Gives dignity when taking patients to and from the shower room.
- Easy to put on by simply putting it over the patient's head.

Shower Cape	
Code	Style
1028169	291
1258919	228



## MOBILITY AIDS

### Aspire Star Patient Mover



Offering a more active transfer solution, users are encouraged to stand up independently, improving circulation, breathing and muscle strength. The Aspire Star Patient Mover is an ideal alternative to a patient lifter for users with appropriate sitting balance, who require transfers over greater distances.

Star Patient Mover	
Code	1461575

### Gutter Frame Rollator Walker



The Days Gutter Walker with Handbrakes has a NEW IMPROVED DESIGN with thicker tubing and larger wheels. Designed to keep the user in an upright position for gait training. Folds for storage. Forearm and hand supports are fully adjustable.

- Height: 1100mm, Width: 630mm
- Weight Capacity: 150kg
- Between Wheels: 580mm

Gutter Frame Walker	
Code	1304534

## Seat Walkers

- Basket and handbrakes
- Adjustable height 830mm – 930mm
- 150kg capacity



Wheelie Walkers	
Code	Colour
1250017	Black
1396781	Red
1396773	Blue

## Days Link Wheelchair



- Self Propelled
- Lightweight, adaptable
- Practical

Days Link	
Code	Spec
1292986	20" 150kg
1284290	18" 130kg

## Days Whirl Wheelchair



- Robust steel frame
- Removable armrests
- Swing away and removable leg rests
- Puncture proof solid tyres
- 18" wide seat with 120kg user capacity

Days Whirl	
Code	Spec
1279130	Self Propelled
1279149	Transit Attendant Propelled

## Days Swift Wheelchair



- Robust steel frame
- Good manoeuvrability
- Fixed Arms Swing away and removable leg rests
- Puncture proof solid tyres
- 18" wide seat with 115kg user capacity
- Ultralight 7.8kg

Days Swift	
Code	Spec
1324810	Self Propelled
1331957	Transit Attendant Propelled

## Sentra Bariatric Wheelchair



- Sturdy, carbon steel frame
- Simple folding mechanism
- Swing in/out leg rests
- Swing back removable armrests
- 200kg max user weight

Days Link	
Code	Spec
1274420	24" 600mm
1273183	22" 550mm



## ID Soft Tilt-In-Space Reclining Wheelchair

- Sturdy steel frame providing extra support and stability in any position for added security
- Infinite position adjustment of the backrest and seat
- Larger wheels for a smoother ride
- 130kg Maximum user weight

ID Soft Tilt-in-Space	
Code	Spec
1297562	400mm 18"

## Footplate Covers

- Sheepskin and velour footplate covers
- Provides protection for shins

Code: 1027635





## WARMTH & COMFORT


### Beurer Electric Heat Warmers



- Cosy and soft – cover made from microfibre
- Overheating protection
- Automatic-off function after approx. 90 minutes
- 3 temperature settings
- Machine-washable at 30°C

Heat Warmer	
Code	Spec
1414909	44x33cm

### Sheepskin Underblanket



- Natural fibres for better temperature balance
- Withstands high washing temperatures
- Urine resistant

Underblanket	
Code	Spec
1172672	30mm

## MEALTIME AIDS

### Felt Clothing Protectors



- Dignified protection for personal clothing
- Highly absorbent with impermeable barrier backing
- Adjustable snap closure
- Full (90 x 45cm) or half length (45 x 45cm) options available

1221490	Half Length Azure Blue
1289659	Full Length Blue Plaid
1287028	Full Length Maroon Plaid



### Rolyan Manoy Built Up Plates

- Plates with sloping bottom are ideal for people who have the use of only one hand
- Dishwasher safe
- Not microwave safe

1053368	Small Plate 220 x 180 mm	Each
1053359	Large Plate 280 x 200 mm	Each
1028347	Plate Sorround	Each



### White Terry Towelling Clothing Protector

1215368	Pop Over Style with Yellow Binding
---------	------------------------------------



### Manoy Economy Scoop Dish

- Moulded low in front and high in back, this unbreakable, round scoop dish is especially useful.
- Non-slip, rubber-padded bottom provides control when using the scoop feature. \* Diameter: 200mm.
- Washes in top shelf of domestic dishwashers.
- Not recommended for microwave use.

1172799	Economy Scoop Dish	Each
---------	--------------------	------



### Rolyan Supagrip Cutlery Range

1184889	Knife	1184904	Tablespoon
1184898	Fork	1395319	Teaspoon

## T.E.N.S THERAPY



### Omron® Deluxe HV-F128

The Omron HV-F128 is the deluxe T.E.N.S device from Omron. With twelve automatic therapy modes the Omron HV-F128 can easily be set to offer either pain relief or muscle stimulation at the touch of a button.

- Provides relief of pain
- Suitable for management of muscular injuries
- High frequency stimulation to relieve pain quickly
- Single channel device

1242006	HV-F128 T.E.N.S Machine
---------	-------------------------



### Omron® Standard HV-F127

Compact in its design the Omron HV-F127 is a simple to use T.E.N.S device. With one touch operation you can easily select one of the nine automatic therapy modes.

- Provides relief of pain
- Suitable for management of muscular injuries
- High frequency stimulation to relieve pain quickly
- Single channel device

1241591	HV-F127 T.E.N.S Machine
---------	-------------------------



### Omron® T.E.N.S Electrodes

- Long Life Pads
- No messy gels required
- Suitable for use with Omron HV-F128 & HVF127

1170373	Omron Electrode Pads	Pack 2
---------	----------------------	--------

\* Not to be used in conjunction with electronic implants such as pacemakers. Read Manufacturers instructions before using these products



### ProTens Dual Channel T.E.N.S Machine Kit

Kit contains: one dual-channel T.E.N.S stimulator, carrying case, lead wires, operating manual, 9V batteries, pack of four electrodes 50mm

1280813	Pro T.E.N.S Machine
---------	---------------------



### Pals T.E.N.S Electrodes

To suit Protensens machine

1326163	Valutrode 5 x 5cm Square	Pack 4
1252316	Platinum 5 x 5cm Square	Pack 4
1054365	Platinum 5cm Round	Pack 4
1064058	Platinum 5 x 9cm Rectangle	Pack 4

## NURSING AIDS

### Pelican™ Smoking Aprons

- Designed to protect the smoker from burns if they drop a cigarette
- Made from kevlar and nomex (used by firefighters and racing car drivers)
- Can be washed without losing it's fire retarding properties

Smoking Aprons	
Code	Spec
1258919	Without Side Flaps
1287354	With Side Flaps



## CRASH MATS



### PremiumLift Crash Mat

- Helps reduce impact forces by up to 90%
- Helps provide full coverage of the bed exit area
- 3cm thick foam creates low profile edge
- Tri-fold design for easy storage
- TGA Approved
- Wipe clean vinyl cover
- Dimensions: 1780 x 965 x 30mm (LxBxH)

1453726	1780 x 965 x 30mm
---------	-------------------



### HealthSaver Bedside Crash Mat (Droppies)

- Made from high-impact foam with an easy-to-clean vinyl cover.
- Non-slip material ensures that the floor mat remains firmly in place
- Sloping edges make it easy to manoeuvre a wheelchair over the fall mat.
- Dimensions: 1800 x 1000 x 50mm (LxBxH)

1325485	1820 x 1000 x 50mm
---------	--------------------



## SLIDE SHEETS



### Silicone Turning Slide Sheets

- These special slippery sheets are used for turning residents in bed and providing easier manual handling
- Made from washable, drip-dry, light-weight fabric that is silicone coated

Slide Sheets		
1240413	Roll Blue 50m x 1.5m	Each
1274037	Double Red 2 x 1.5m	Each
1396307	Double Green 2 x 1.5m	Each



### Slide Sheet Holder

1362720	Blue	Each
1334638	Red	Each
1293311	Colour Varies	Each



### Wendy Lett Slide Sheets

The WendyLett base sheet is a unique sliding sheet that appears to be a white and grey sheet with narrow white stripes, and is available in many different sizes. The sliding takes place sideways along the direction of the stripes.

1331353	4 Way Draw King Single 200x200cm
1393650	4 Way Draw Single 140x200cm
1331396	Base Sheet King Single 105x200cm
1331361	Base Sheet Single 90x200cm

## PATIENT SLIDE

### Pat-slide® All Purpose Resident Slide



- An efficient, cost effective and comfortable method for general resident transfers
- Safe and easy to use
- Static resistive
- Light weight, radiolucent, easy-to-clean and maintenance free

1078788	Patient Slide 1525mm (L) x 635mm (W)	Each
---------	--------------------------------------	------

## WALK BELTS

### Pelican™ 4 Handle Walk Belts



- Four handle and rubber back make it safe, supportive and easy to use
- Useful when transferring residents

1154351	Small Red	Each
1154360	Medium Yellow	Each
1154342	Large Green	Each
1150183	X-Large Blue	Each

## SIT & SLIDE PAD Pelican™ Sit & Slide Chair Pad



- Allows resident to sit down at the front of a chair and then slide hips back into the chair
- One or two nurses can give assistance by pulling the pad back with the lifting loops
- The resident can also be helped to stand by carer pulling the webbing handles

1125841	Slide Pad for Chairs with Handles	Each
---------	-----------------------------------	------

## RESIDENT TRANSFER BELT

### Pelican™ Soft Transfer Belt



- Four large padded vertical handles
- Easily adjusted for size to suit a variety of residents
- The soft padded belt is slip resistant

1164358	Standard 80-160cm	Each
---------	-------------------	------



## EMERGENCY RESIDENT TRANSFER

### Pelican™ Patient Rescue Sheet

- Kept in position under the mattress with strong elastic corner straps
- Daily bed making is as normal
- In case of fire or emergency, the webbing straps are pulled out and quickly done up around the resident
- The mattress and the resident can then lowered onto the floor for very fast and safe evacuation, without staff having to carry them
- 300kg Capacity

1027949	Rescue Sheet for Vinyl Only	Each
1249740	Rescue Sheet for Vinyl & Carpet	Each

## HOIST SLINGS

### Pelican™ Dex Full Body Sling with Loops

- General purpose designed sling with head support
- All webbing is sewn on the outside of the sling to prevent pressure


All Day Green Nylon	
1137403	Large
Blue Polyester Mesh	
1320017	Small
1209164	Medium
White Nylon	
1416634	Small Red
1416642	Medium Yellow
1137403	Large Green



### Pelican™ Premier Lifting Slings


- Full body sling with soft padded split legs and padded head support
- Designed to give a variety of lifting positions to enable residents to be lifted in a reclined mode or in a sitting upright position

White Nylon	
1027840	Medium
1102785	Large
1281577	X-Large
1155789	Adult Head & Shoulder Small
All Day Green Nylon	
1376802	Small
1366459	Medium
1366467	Large
Blue Polyester	
1279084	Small
1279092	Medium
1279106	Large




### AidaCare Aspire Deluxe Slings


- Australian standards tested and registered with the TGA
- Intuitively designed, this sling offers generous underarm padding and a reinforced wide-cut back for comfort
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting



General Purpose



Hygiene Access



Padded Standing

Item #	Aspire Deluxe Swings	Size	Load
1415409	General Purpose Sling with Head Support	Small	200kg
1425587	General Purpose Sling with Head Support	Medium	200kg
1425609	General Purpose Sling with Head Support	Large	300kg
1425595	General Purpose Sling with Head Support	X Large	300kg
1415220	Hygiene Access Sling with Head Support	Small	200kg
1415204	Hygiene Access Sling with Head Support	Medium	200kg
1415212	Hygiene Access Sling with Head Support	Large	300kg
1415239	Hygiene Access Sling with Head Support	X Large	300kg
1425560	Padded Standing Sling	Medium	200kg
1425552	Padded Standing Sling	Large	300kg

## HOISTS



### Alphacare Lopital Reflex Electric Shower Commode

**\*See page 5 for more information**

Code	Max Weight	Seat Height	Seat Width	Overall Width
1435310	135kg	490-1001mm	480mm	660mm

### Aspire 200 Steel Lifter

- Designed to meet the rigorous demands of professional and facility environments
- Ergonomic handles provide multiple grip points to suit carers of different heights and strengths
- Powered leg function opens and closes legs for lifting stability or accessibility
- Strong powder coated steel frame with industry leading Linak electronics
- Critical safety features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- Warranty - Frame 5 Years, Electronics 3 Years



Item Code	Description	Max User Weight	Leg Gait	Lift Range
1373056	Aspire 200 Lifter	200 Kg	580 - 900mm	535 - 1725mm

### Aspire 200 Steel Standing Lifter

- Bullhorn Lifting Handles for patient stability and 6 sling attachment points
- Ergonomic Handles Provide multiple grip points to suit carers of different heights and strengths
- Ergonomic Handset The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Powered Leg Function Opens and closes the legs for lifting stability or accessibility
- Warranty - Frame 5 Years, Electronics 3 Years, Battery & Handset 1 Year



Item Code	Description	Max User Weight	Leg Gait	Lift Range
1391054	Aspire 200S Standing Lifter	200 Kg	575 - 936 mm	850 - 1740 mm

### Aspire A205 Aluminium Patient Lifter

- The Aidacare Aspire A205 is a compact and versatile aluminium mobile lifter with powered legs that is ideal for home, community and facility environments.
- Transportation or Storage: Is quick and convenient with tool free disassembly
- Ergonomic Handset: The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Adjustable Height: Three adjustable height settings provides the caregiver with the flexibility to adjust the height range to suit the specific lift and patient scenario



Item Code	Description	Max User Weight	Leg Gait	Lift Range
1422766	Aspire A205 Lifter	205 KG	547 - 1015mm	395 - 1805 mm



### Aspire Inline Weigh Scale

- Designed to be fitted on Aidacare Aspire Patient Lifters
- Allows for an accurate weight measurement to be taken by a carer when a patient is suspended in a sling
- Large weight capacity of up to 320kg
- Easy to read LCD display

Item Code	Max User Weight	Warranty
1415727	320 kg	2 Years



## Alphacare Smart Stand Patient Transfer Aid

- Foot Operated Leg Spread
- Swing Away Seat
- Waterproof Padded Seat Pads
- Rear Braking Castors
- Optional Pelvic Belt
- Durable Powder Coated Frame
- 182kg Max User Weight

Code	Max Weight
1435329	182kg

## Aspire Swift Patient Mover

- Four user contact points (foot platform, knee pad, seat pads and grip handle) promote stability, support and comfort for the user - resulting in greater confidence in the transfer
- Rigid foot platform - provides a stable surface for the user to stand on and support their own weight

Item Code	Description	Max User Weight	Total Width	Total Height
1396528	Switch Patient Mover	185kg	560-710mm	1060mm



## SENSOR BEAM

The ultimate early warning for fall preventions.

- Utilises motion detection
- Creates an invisible beam above or beside a bed/chair or across a doorway
- Generates a nurse call alarm when beam is interrupted
- Connects to existing nurse call system
- Great for residents who have become aware of floor mats
- Will not false alarm
- Transferable from room to room

1281321	Nurse Call Alarm Sensor Beam w 3 Metre Call Bell Pole & Base
1369121	Nurse Call Alarm Sensor Beam Call Bell with Wall Mount Brkt
1373153	Nurse Call Alarm Sensor Beam with Small Pin RCA
1402005	Sensor Beam w Phone Style Call Plug,Bell,Pendant,Pole & Base




# PRICE MATCH POLICY

At SSS Australia we regularly review our prices to ensure we remain competitive, and you always get a great deal on your medical supplies.

If you see a competitor quote a lower price on the same stocked item, please give us a call and take advantage of our price match policy.





## NURSE CALL BED ALARMS


Out of bed alarms are an inexpensive and highly effective method of reducing falls by monitoring the out of bed movements of at risk residents.

The Out of Bed Sensor Mat can be easily fitted to your existing nurse call system. The Sensor Mat interface control replaces your existing nurse call pendant (nurse call push button on a cord) with a new pendant that also connects the Out of Bed sensor mat to your nurse call system. The new pendant & interface control feature an on/off switch for your convenience.

Sensor Mats	
1415506	Stand Alone Out of Bed Sensor Mat Grey 900 x 600mm
1335413	Stand Alone Out of Bed Sensor Mat Grey 1200x600mm
	1415506 and 1335413 Requires Interface Control



Interface Controls for Out of Bed Sensor Mats	
1281305	Interface Control & Pendant A Stereo (2 Rings)
1332678	Interface Control & Pendant A Mono (1 Ring) No Nurse Call
1321560	Interface Telephone Plug
1281410	Interface Control & Pendant B 1 Pin



Please speak to your personal account manager about connection requirements for your facility.

## SAFETY SOCKS

### Ankle Height Safety Socks

- Non-slip grip spots on top and bottom of socks
- Bright orange for easy identification of residents at risk

1320157	SafeStep Safety Socks Regular Sml-Med (Sizes 2-6) Orange
1320173	SafeStep Safety Socks Large (Sizes 7-11) Orange
1320181	SafeStep Safety Socks Extra Large (Sizes 12-15) Orange
1323040	Safety Socks Med (Med-Lge Sizes 7-11) Black Gripperz
1323059	Safety Socks Small (Sml-Med Sizes 5-7) Black Gripperz



### Foot Slipper Safety Socks

- Long. Double sided grip.

1331175	Foot Slipper Safety Socks Long Double Side DkGreen Bariatric
1331183	Foot Slipper Safety Socks Long Double Sided Red Sml/Med
1331191	Foot Slipper Safety Socks Long Double Sided Red Lge/XL
1331205	Foot Slipper Safety Socks Long Double Sided Dk Green Sml/Med
1331213	Foot Slipper Safety Socks Long Double Sided Dk Green Lge/XL



## FRIDGES



### Avem Quirks MEDISAFE PLUS Purpose built QCPP accredited vaccine fridges

- Specifically designed for vaccine storage with clear glass doors and temperature controls
- Reliably maintains the temperature between 2°C & 8°C
- 24 hour product support
- 2 year warranty

**INCLUDES 1 YEAR OF CLOUD BASED MONITORING & AVEM QUIRKS BOX**

1279211	Medisafe Gen2 Glass Door Vacc Fridge 80 Ltr w SIMCard Cloud
1244486	Medisafe 311Gen2 Glass Door Vacc Fridge 306L w SIMCard Cloud
1244495	Medisafe Gen2 Glass Door Vacc Fridge 381Ltr w SIMCard Cloud
1245400	Medisafe Spare Shelf for 306ltr & 381ltr Vaccine Fridge
1244501	Medisafe Auto Dialler Fridge Option
1385917	Medisafe Fridge AQ Box Cloud Based Temp Monitoring Device
1329839	Logtag Software & Docking Station Only
1329847	Logtag Temperature Data Logger with Probe

### UPS and Battery Back Up Pack



- UM-IP2000 is specifically designed for Medical refrigeration providing extended power back up during a power outage.
- 3 year warranty

1436813 Medisafe Fridge UPS and Battery Pack



### NATION WIDE FIELD SUPPORT

Put your inspection, calibration and testing requirements into safe hands.


- Programmed maintenance on your vaccine fridge
- Equipment calibration
- Equipment tag and testing
- Dedicated team of over 85 technicians

**For service & support phone:  
1800 658 120**





## PREMIUMLIFT ULTRA LOW BED



## CodaCare

CodaCare is 100% Australian owned. Their story starts way back in 1967, and to this day they continue to build on their heritage of user-centric design, quality local manufacturing, importing and outstanding client service. As leading hospital & aged care residential furniture suppliers, they believe in providing long term value for money for Australia's hospitals, rehabilitation centres, aged care facilities and other healthcare organisations.

### Available in single and king single

The PremiumLift Ultra Floorline Bed is the ultimate in low level beds, designed in Germany to the higher quality specification. It offers a variety of functions including Trendelenburg and cardiac chair position allowing for excellent care to be provided. It is an optimal floor level bed solution for use both for home care and aged care facilities. Its minimum height from the floor is just 6.7cm ensuring complete patient safety, should the user be inclined to roll whilst in bed, and potentially fall, the bed reduces impact and allows the user to safely roll onto the floor, or a fall matt.

1427962	PremiumLift Ultra Low Bed Beech King Sgl 7cm - 64cm SWL 185kg
1427970	PremiumLift Ultra Low Bed Beech Single 7cm - 64cm SWL 185kg

They use a range of materials to produce furniture which includes hospital grade compact laminate, a marine grade board that comes with a 10 year warranty. HPL is impact and scratch resistant and does not fade over time.

When you buy Australian made furniture, not only are you supporting local manufacturers and retailers, you are choosing furniture made from quality and sustainable materials which are manufactured under the guise of quality management system such as ISO 9001 and Australian health and safety requirements.



### PremiumLift Premium Mattresses

- Three layer foam construction
- Super soft memory foam top layer
- High density medium firmness center support layer
- High density firm base layer to prevent bottoming out
- Safe working load of up to 140kg - 15mm thick
- Therapeutic working load of up to 150kg
- Medium risk pressure care properties
- Extra supple Two-Way Stretch cover

1454404	PremiumLift Premium Mattress SINGLE (1930cm x 900cm x 150cm)
1453718	PremiumLift Premium Mattress KING SINGLE (2030 x 1060 x 150cm)



## BEDROOM ACCESSORIES

### CodaCare Premiumlift Overbed Table



- Sturdy steel frame powder coated
- Spring adjustable with locking hand wheel
- Moisture resistant laminated top, 800 x 400mm (18mm thick)
- Impact resistant PVC matching edging
- 50mm twin wheeled castors (inclusive of 2x lockable castors)
- Height adjustment range 740 to 1050 (mm)
- Safe Working Load: 25kg evenly distributed
- Australian Made

1427954	Overbed Table Beech C-Base Raised Edge 225T
---------	---



### ThermoForm Over Bed Table

This overbed table features a solid "C" base frame and flat top, ideal for use in the hospital and aged care environments. Supplied with a standard willuna white (beige) laminate top with a vinyl wrap around the outside, or with the deluxe "thermoform" willuna white (beige) top which is heat resistant and ideal for hospitals higher requirement needs.

- Thermo heat resistant top with raised edge
- Under-table strata for improved strength
- Lockable castors

1283324	Willuna White Overbed Table
---------	-----------------------------

### Mobile Stainless Steel I.V. Pole



- Five star solid steel base
- Finished in hammertone gunmetal grey powder coat with black 50mm twin wheel castors
- Chrome plated stem with key wheel for height and tension adjustment

1087357	I.V. Stand Stainless Steel Adjustable Height
---------	--

### CodaCare Bedside Locker 3 Draw



- Lockable top drawer providing secure storage for medication and/or valuable personal belongings
- Large, easy-to-grab handles
- Designed with the aged care environment in mind - strong and durable
- Fitted with 4 castors
- 800mm (h) x 450mm (w) x 500mm (d)


1264952	3 Drawer Bedside Locker Beech
---------	-------------------------------

### Bed Rail Protectors



1373617	50x166cm slip over
1027377	50x166cm velcro on waterproof

### Bed Cradle Blanket Support



- Designed to fit under the mattress to keep linen clear of the feet and legs

1096042	Bed Cradle Blanket Support Epoxy
---------	----------------------------------

## MEDICATION TROLLEY

### Mobile Medicine Trolley 4 Drawer White with Accessories

#### Drawer Sizes

- 4 drawers x 110mm
- 1 drawer x 140mm

#### Features

- Non slip mat
- Luxurious push handle
- Soft closing drawers
- Dimensions: 650 x 480 x 900mm
- G.W 42kg
- Extension Table
- Waste Bin
- 304 Stainless Steel Wire Basket
- Central or Seals Lock
- Castors with brakes
- \* All accessories included

2 year Warranty



1337521	Trolley Medicine 4 Drawer White with Accessories
---------	--



1451154	Lavender
1453076	White
1450735	Light Blue

### Viva Narrow Cart Aluminium 34"

The Aluminium Narrow Cart has many of the same features as the traditional Viva Narrow cart, yet with sleek lines, contemporary handles and a lightweight body for a greater ease of use.

- ABS top with slide out work surface
- Central key lock control
- Push handle
- 4 x 100mm twin swivel castors; two locking
- Soft close drawers
- Colour: Lavender
- Width: 510mm
- Depth: 542mm
- Height: 977mm
- 2 Draws: 76mm
- 2 Draws: 156mm
- 1 Draw: 234mm
- Weight: 49kg



1453254	Viva GC1040 3 Drawer 690x520x930
1452959	Viva GC1050 3 Drawer 690x520x1010
1453203	Viva GC1060 4 Drawer 690x520x1085
1452967	Viva Accessory Pk for Isolation Cart

### Viva Isolation Cart Yellow GC1050 3 Drawer

The bright yellow colour of the Viva Isolation Cart is synonymous with the high risk infection control areas in the facility. Safe and dust free storage of clean gowns, caps, gloves etc provide peace of mind for staff and visitors.

- ABS top with slide-out work surface
- Divider in two drawers
- Central key lock control
- Push handle
- 4 x 125mm Twin swivel castors; two locking, one directional
- Soft close drawers
- Colour: Yellow
- Width: 690mm, Depth: 520mm, Height: 1010mm
- 3 Draws: 234mm
- Weight: 55kg

Available in two additional sizes with 3 or 4 drawer options.

Also Available - GZ2680 - Isolation accessory package

- Waste container set and lid
- Multi-storage basket
- Glove Dispenser



## STAINLESS STEEL TROLLEYS

### Econo Trolleys

- Welded stainless steel construction with safety fold to edges
- High grade 304 stainless steel
- Unique self-closing drawer runners
- Smooth 95mm noise-free castors with front brakes
- Draws 125mm High
- Magnetised drawers to prevent opening



1317458 Trolley S/Steel 2 Drawer Shelf & Rails  
50 x 50 x 90cm Econo



1332465 Trolley S/Steel Shelf & Rails/No Drawer  
50 x 50 x 90cm Econo



1336096 Trolley S/S 2 Drawer Side by Side & Shelf  
80x50x90cm Econo



1335464 Trolley S/Steel Single Drawer 50 x 50 x 97cm with Lock 1451



1309609 Trolley S/Steel 1 Drawer Shelf & Rails  
50 x 50 x 90cm Econo



1332783 Trolley S/Steel 4 Drawer/ Shelf & Rails  
50 x 50 x 90cm Econo



1323024 Trolley S/Steel 3 Drawer/ Shelf & Rails  
50 x 50 x 90cm Econo



### OXYGEN TROLLEYS

1093194 Suits size C & D Cylinders Retractable  
1228043 Suits size E Cylinder



### Stool Gas Lift Grey with Foot Rest

- Adjustable height: 49cm - 63cm
- Cushion diameter: 35cm
- Foam thickness: 6cm
- Base: 56cm

1113591 Stool

## FURNITURE FOR THE WHOLE FACILITY

SSS Australia is able to assist with all furniture needs for your facility — from bedroom to lounge room.

Call our friendly Customer Support Team on **1800 777 518** for more information about our furniture range.

## COMMON AREA FURNITURE

### Aspire Mobile Air Chair

- Ergonomic Attendant Grip Handles
- Swing Away Backrest Wings & Armrests. Allows for easy patient side transfers
- Designed with the aged care environment in mind.
- Lockable Castors with Directional Control
- Unique design promotes skin immersion and envelopment, redistributing skin pressure, reducing the risk of pressure injury associated with prolonged sitting or lying

Max User Weight: 180kg

1368516 Aspire Mobile Air Chair Large



#### Footrest for Aspire Mobile Air Chair

- Comfortably padded footrest
- Additional add-on for Aspire Mobile Air Chair

1368524 Footrest Padded for Aspire Mobile Air Chair



#### Meal Tray for Aspire Mobile Air Chair

- Perfect for meal times
- Additional add-on for Aspire Mobile Air Chair



1368532 Meal Tray for Aspire Mobile Air Chair



### Care Quip X6 Deluxe Mobile Air Bed/Chair W- Tray & Footrest

- Ideal for dependent clients who are confined to a chair for long periods
  - Pressure relieving memory foam, contoured backrest with lumbar and lateral support, angled armrests and adjustable wings and neck pillow ensure maximum comfort and postural support
  - Integrated air cell in seat for premium pressure relief
  - Tilt in space recline action and high leg rest elevation for reduced fluid retention
  - Designed for use with standing machines and transfer hoists
  - Large push handle and simple recline mechanism makes movement and operation safe and easy for the carer
  - Resilient European four-way stretch vinyl is fluid resistant, fire retardant and vapour permeable with antibacterial and non-allergenic properties
- Seat Width: 500mm
  - Back Height: 650mm
  - Seat Depth: 500mm
  - Overall Width: 750mm
  - Arm Height: 180mm
  - Max User Weight: 180kg

1451669 Care Quip X6 Deluxe Mobile Air Bed/Chair W- Tray & Footrest



#### Aspire Adjustable Day Chair

- High back orthopaedic chair that offers adjustability, comfort and support
- 160kg safe working load
- Latte vinyl colour
- Depth adjustable
- Contoured backrest

1393057 Adjustable Day Chair





#### Aspire Waterfall Day Chair

- High back orthopaedic chair that offers adjustability, comfort and support
- 200kg safe working load
- Ink vinyl colour
- Padded seats and armrests
- Height adjustable
- Contoured waterfall backrest - segmented and adjustable backrest cushion provides comfort and postural support

1391402 Waterfall day chair

## MATTRESSES



### Icon Pressure Care Mattress

- Profile stabilising sides enables ease of care
- Triple layered, premium pressure care core designed to provide optimal pressure care whilst providing support and comfort
- Premium visco grade, pressure sensitive memory upper minimises resident interface pressure by maximising contact surface area
- Water-based evaporative lamination - completely odour free

#### Optional Extras Include:

- Safety sides - Soft safety sides provide unobtrusive support to help prevent high risk falls residents rolling out of bed
- High frequency sealed and welded cover seams - Provides maximum infection control and facilitates cover cleaning and care

#### Icon I3100

Code	Length	Width	Depth	Notes
1365053	1900mm	880mm	150mm	
1370103	1960mm	880mm	150mm	
1348485	1980mm	880mm	150mm	
1372467	1980mm	1050mm	150mm	*suits user up to 280kg King Single
1358103	2030mm	1050mm	150mm	*King Single
1359460	1940mm	890mm	150mm	*w 245mm concave

#### Icon IM1000 Bariatric

Code	Length	Width	Depth	Notes
1404555	1980mm	1050mm	180mm	*suit users up to 350kg
1444158	1980mm	880mm	150mm	*120kg-280kg

### Icon Pressure Care Mattress Budget

- Triple foam layer core with pressure sensitive visco memory upper for immersion and pressure care
- 4-way Stretch, Vapour Permeable European Pressure Care Upper
- Heavy Duty PVC Base

#### Icon I21000 Budget

Code	Length	Width	Depth	Notes
1364154	1900mm	880mm	150mm	
1358758	1980mm	880mm	150mm	

### Sovereign Superior Pressure Mattress

- Triple layered core - premium medical grade, heavy weight foams
- Heel gradient - minimises pressure loading on critical heel areas
- Ventraflow inner core technology - internal air circulation provides over 3 times greater internal air flows minimising internal heat entrapment and maximum moisture evacuation
- Airo-form visco pressure sensitive memory upper - minimises resident interface pressure by maximising contact surface area

#### Sovereign S41000

Code	Length	Width	Depth	Notes
1380958	1980mm	880mm	150mm	
1394665	2030mm	1050mm	150mm	King Single





The Forté Assure Evacuation System eliminates the annoyance of traditional mattress evacuation systems through its totally integrated unique design. With all evacuation straps self-contained within inside the mattress cover, there is no safety or damage issue of straps falling loose and becoming entangled with nearby medical equipment or care staff attending to the patient. With all seams High Frequency Sealed and Welded, the Assure Evacuation System does not compromise on hygiene and infection control. Staff safety and risk of injury is mitigated as the Assure System significantly reduces force required to evacuate a patient across a range of floor surfaces. Unique foam core modifications allow the patient to be encapsulated and further protected against harm.



## Evacuation Pressure Care Mattress

- Custom Foam Core Modification
- Allows foams to encapsulate patient for protection
- Low Resistance Base Sheet
- Tapered to allow grip to bed base while in normal use
- Significantly reduces friction forces between mattress and varied floor surfaces to make evacuation procedure easier
- External Seams High Frequency Sealed, Welded and Fused for added infection control
- Prevents Ingress of Fluid through Seams
- Compliant with Infection Control Requirements
- Fully Integrated & Contained Evacuation Straps
- Wipe Clean Pocket with High Frequency Welded Zipper for Infection Control
- Easy to Access, do not need to lift mattress to access
- Pockets Integrated within foam core
- Dual Evacuation Pull-Point
- 4 x Integrated Robust Heavy Duty Carry Handles
- 100% Wipe-Clean, Impervious to Fluids
- Stainless Steel, Non-Corrosive Nyloc Attachment

Codes	Length	Width	Depth	Notes
1454064	1980mm	880mm	150mm	Medium Risk
1454129	1980mm	880mm	150mm	High Risk



## Zephair Pressure & Palliative Care Mattress System

High Level Acute Pressure Care and Comfort Mattress System  
The ZEPHAIR Mattress range provides true leadership in powered patient pressure care. Tried and trusted in facilities from ICU Wards to Palliative Hospice use, this pressure care system will provide the answer to all acute care needs, whilst providing unparalleled comfort and staff Primary Care compliance.

Key clinical attributes include ventilated immersion upper layers, deep 24 cell alternating integrated inner mattress core, firm base foam layer, quiet efficient air flow generator with easy set weight calibration and dual stretch waterproof breathable cover.

The ZEPHAIR Pressure Care Mattress range includes models suited for standard hospital bed sizing, wide & bariatric requirements. The quality and features of the Zephair range allow confident usage in the Hospital, Aged Care and Spinal environments. Select the Zephair Pressure Care System and achieve unsurpassed reliability and consistency of care.


- 24 Cell, Alternating Air System
- 125mm depth facilitates immersion
  - Ultra-flexible, soft-touch material
  - 24 Cell System minimises patient disturbance
  - Static head support region



Zephair ZA1000				
Code	Length	Width	Depth	Notes
1395815	2000mm	900mm	220mm	Single
1395823	2000mm	900mm	220mm	King Single

**5 year foam core warranty**  
**2 year air cell warranty**  
**1 year pump warranty**

## MATTRESS WITH PRESSURE REDUCING INSERT



### Mattress with Customised Roho® Insert

- Designed as a high level product for pressure reduction
- Suitable for use with High Risk and Palliative Care
- Available with any Forte static support surface
- Custom made to resident/facility requirements

Speak with your Personal Account Manager about options to suit your needs

24 month warranty

Code	Length	Width	Depth	Notes
1396161				Single Mattress
1297708	860mm	51mm		Roho Insert

## ALTERNATING MATTRESSES



Shown without multi-stretch cover

### Smik Care™ Alternating Overlay Pressure Mattress System

- Assists carers to provide maximum pressure relief
- Multi stretch vapour permeable cover
- Equipped with alternating, static and auto firm (nursing) modes
- Auto firm (nursing) mode maximises the internal pressure
- Retainer air cell design prevents over inflation
- Max user weight: 150kg

Code	Length	Width	Depth	Notes
1365657	1950mm	950mm	130mm	Single Mattress with pump
1374303	1950mm	1070mm	130mm	King Single with Pump



### Forte Airmonte A5000 Alternating Air Mattress Overlay

- For medium to high risk
- Two-way stretch, vapour permeable cover
- 160kg capacity

Code	Length	Width	Depth	Notes
1402862	1960mm	860mm	125mm	150kg capacity



### Forte Airmonte A8000 Full Replacement Air Mattress System

- High quality pressure redistribution surface suitable for residents with medium to high risk of Pressure Injury
- Specially designed cells incorporated with a low pressure high flow pumping system, enables maximum pressure redistribution
- 180kg capacity

Code	Length	Width	Depth	Notes
1402854	1980mm	870mm	150mm	200kg capacity
1426257	2010mm	1080mm	150mm	King Single



## MATTRESS UNDERLAY


### Alternating Mattress Underlay

Code	Length	Width	Depth	Notes
1271948	1980mm	880mm	50mm	Single
1376942	1980mm	1050mm	100mm	King Single

## PRESSURE CUSHIONS

### EquaGel® Cushions

- EquaGel is a dry-polymer gel, it has unbelievable tensile and compression strength and incredible durability.
- 2 Layer Gel technology- top layer is a thin wall with tight grid for comfort, the bottom layer is thicker for deep posture support
- Grid design spreads & reduces pressure over a wider area
- Contoured pelvic area of the cushion for additional equalizing pressure distribution
- Black cushion cover is fluid resistant and vapour permeable
- Moderate General Cushion - 2" tall maintaining 1 1/4" gel below pressure points
- Advanced Protector Cushion - 2.5" tall maintaining 1 3/4" gel below pressure points



Moderate General Cushion		
Code	Length	Width
1327712	16"	18"
1369709	20"	16"
1334360	18"	18"

Advanced Protector Cushion		
Code	Length	Width
1327720	16"	18"
1334379	18"	18"
1418203	18"	16"
1418157	16"	16"

Adj Advanced Protector Cushion		
Code	Length	Width
1369717	18"	16"
1381601	22"	18"
1426354	24"	20"

### Roho® Quadro Alternating Cushion

- ISOFLO® Memory Control offers shape fitting capabilities while the user is seated
- Quick and easy on-demand adjustment to maximize function
- Facilitates blood flow by fitting, matching and tracking the shape of a person
- Cell size: 9 x 9cm

Code	Desc.	Length	Width	Depth
1286048	HP Cushion	410mm	410mm	90mm
1328018	Cushion cover for Roho Quadro			



### Ultimate Bliss Premium Cushion

- A combination of high density and memory foam helps provide pressure reduction and comfort for the at risk resident
- Cover features a vapour permeable and two-way stretch waterproof top with a non-slip base, entire cover can be cleaned for easy maintenance

Code	Desc.	Length	Width	Depth
1331566	Cushion	450mm	450mm	30mm
1272201	Cushion	450mm	450mm	90mm
1272210	Wheelchair Cushion	400mm	400mm	70mm



### Treat-Lite® Air Cushion

- An air cushion without a pump - no adjustments, no maintenance
- High risk cushion which can have the inner layers removed for varying depths
- Can be machine washed

Code	Desc.	Length	Width	Depth
1387553	Cushion	450mm	450mm	30mm
1387561	Cushion	500mm	450mm	90mm



### Forte Sonnet 2 Core Cushion

- High performance medical grade foams to provide sustained patient care and comfort

Code	Desc.	Length	Width	Depth
1418939	Cushion	460mm	400mm	70mm
1418920	Cushion	460mm	460mm	70mm
1419455	Pressure relieving heel cushion			






## BODY PROTECTION

### DermaSaver®

- MicroSpring textile for great protection
- Reduces friction and pressure against the skin
- Many designs available – ensures an ideal product to suit your residents' needs



Finger Separator		
1325493	One Size Fits All	Each
Forearm Tube		
1388088	Small 18-25cm	Each
1388096	Medium 25-30cm	Each
1388118	Large 30-38cm	Each
1388126	Extra Large 38-46cm	Each
Full Leg Tube		
1325752	Extra Small 18-23cm	Each
1325744	Small 23-30cm	Each
1325736	Medium 30-38cm	Each
1325728	Large 38-48cm	Each
1325701	Extra Large 48-58cm	Each
Heel Protector Stay Put		
1325531	Small 28-30cm	Each
1325523	Medium 30-37cm	Each
1325515	Large 37-39cm	Each
1325507	Extra Large Over 39cm	Each
Heel Protector with Toe Cover		
1376608	Small 28-30cm	Each
1387952	Medium 30-35cm	Each
1387960	Large 35-39cm	Each
1387979	Extra Large Over 39cm	Each
Ankle Bumper		
1377159	Medium 23-28cm	Each
Pillow Case		
1388134	Pressure Reduction Case 66x51	Each

Knuckle Protector		
1388169	Small/Medium Below 20cm	Each
1388177	Large/Extra Large Above 20cm	Each
Oxygen/Feeding Tube Cover		
1325760	Tube Covers	Pk 10
Palm Pillow		
1388142	Small	Each
1388150	Large	Each
Shin Tube		
1387901	Extra Small 18-23cm	Each
1387928	Small 23-30cm	Each
1336452	Medium 30-38cm	Each
1387936	Large 38-48cm	Each
1387944	Extra Large 48-58cm	Each
Toe Tubes		
1388185	10 Pack (3 sml / 5 med / 2 lge)	Pack



### CareProtect Pedi Boot

- A unique device to distribute pressure for residents with a pressure injury risk
- High density foam provides low load-bearing effect
- Open toe allows foot to breathe
- Seams on exterior for less irritation
- Reusable and washable

Code	Length	Notes
1319671	Size 1: Green 27-31cm	Pair
1319698	Size 2: Blue 32-36cm	Pair
1319701	Size 3: Yellow 37-42cm	Pair



sutherlandmedical





## BODY PROTECTION

### Gelbo Elbow / Heel Protector

- Mineral oil-based gel formula: softens and moisturises dry cracked skin
- Reduces pressure, friction and shear forces

1240306	Medium	Pair
1240315	Large	Pair
1240324	Extra Large	Pair



Machine Washable



### Sheepskin Heel and Elbow Protector

- Hi-temp washable
- Adjustable velcro fitting

1184190	Green One Size Fits Most	Pair
---------	--------------------------	------



### Care Quip Silicone Fiber Heel Protector

- Padded heel cushion, provides comfort and assists in the prevention of pressure build-up in the heel area

1180772	Blue One Size Fits Most	Each
---------	-------------------------	------

### Protect-A-Limb®

- Prevents skin tears and damages
- Protection against friction to vulnerable limbs

1268989	Small	Pair
1268998	Medium	Pair
1269019	Large	Pair



### Comfy Ears®

- Cushioned protection for Oxygen Cannula Tubing
- Can be trimmed to required size

1226209	Large Protectors	Each
---------	------------------	------



## HEEL WEDGES

### Pressure Relieving Heel Wedge

Code	Size	Notes
1225640	43x22x9cm	Single
1216033	86x46x9cm	Double
1215878	45x18x9cm	Single Concave



## HEEL ELEVATOR

### DermaSaver® Heel Elevator



- Designed to elevate the heel and cover and protect both of the bony prominences on either side of the ankle

Code	Size
1325949	Heel Elevator (one size fits most)



## HEEL PROTECTOR WITH WEDGE

### Heelmedix Standard Heel Protector

- Completely off-loads the heel and helps redistribute pressure
- Open heel design allows for airflow and easy monitoring
- Straps pull down on sides to help protect against foot drop



1375814	Standard Heel Protector with Wedge
---------	------------------------------------


## BODY WEDGES

### Pressure Relieving Body Wedge

1217729	Wedge 60x30cm angled height 9cm to 14.5cm
---------	---



## SHEETS & PILLOWCASES



### Fitted Sheets

- Designed to encase the mattress and overlay
- Allows the body to conform to the pressure reduction features of the support surface - without hammocking

1263721	PolyCotton Stretch Knit Fitted Sheet 90 x 200cm	Each
---------	---	------

### Excel Flat Sheets & Pillowcases

- Made to Australian standards
- Will not pill
- Only 2% shrinkage

1274912	50% polyester 50% cotton pillowcase 78x50x14cm	Box 100
1274903	50% polyester 50% cotton single sheet 180x290cm	Box 25

## BLANKETS

### Ultra Plush Microfibre Antished

1419544	Charcoal 300gsm	Each
1419552	Mushroom 300gsm	Each

### Celleweave Blanket



100% Cotton Preshrunk White

1337025	Single 180x230cm	Each
---------	------------------	------

### Coral Fleece Large Single

1391410	Walnut 180x245cm	Each
1391437	Chamios 180x245cm	Each

## BATH TOWELS & FACE WASHERS



### Pearl Bathing Range

- Designed to withstand commercial laundries and perform for an extended life
- 100% cotton carded pile, soft to the touch and highly absorbent

Bath Towels		
1308890	155x71cm White 550gsm	Ctn48
Face Washers		
1291092	33x35cm White 56gsm	Pk 50
Hand Towel		
1338803	40x60cm White 120 gsm	Each

### Haines Premium Disposable Bath Towel

1332066	150x80cm White	Pk 10
---------	----------------	-------

### Other Bath Towels & Face Washers

Bath Towels		
1338781	68x137cm White 470gsm	Ctn32
Face Washers		
1290908	28x28cm White Small	Each
1338811	33x33cm White 55gsm	Each
1316443	33x33cm Navy 50gsm	Ctn100

## PILLOWS



### Healthcare Pillow with Smart Barrier Cover

1326953	Pillow waterproof 45x70cm	Each
---------	---------------------------	------



### Vinyl Pillow Protector

1029522	Plastic Zippered	Each
---------	------------------	------

## LINEN PROTECTION



### Pinkie Pads

- Replaces macintoshes and drawsheets
- Large size Pinkie Pads can absorb up to 1L of fluid – with or without tuck-ins
- Small size is perfect to be used as a chair pad or a pillow protector

1288253	Small waterproof 60x43cm (chair)	Each
1161721	Large waterproof 90x85cm	Each
1287346	Large waterproof with wings 90x85cm	Each



### Macintoshes & Drawsheets

- Heavy duty and high quality makes them commercial laundry compatible
- Draw Mac can be used in conjunction with Draw Sheets
- Draw Sheets have a herringbone twill for greater absorbency

Draw Mac		
1218987	Draw Mac White 175x90cm	Each
1289632	Vintex Full Wrap Tuck-in 180x90cm	Each
1275956	Heavy Duty with Wings 100x90cm	Each
Draw Sheet		
1290363	Herringbone Draw Sheet 180x120cm	Each



### Absorbent Pads

Bed Pads		
1286528	2L Capacity with Wings Non Waterproof 100 x 100cm	Each
Chair Pads		
1104189	Small 300ml Capacity 45 x 50cm Tartan	Each
1316877	Large 1L Capacity 50 x 60cm Dark Blue	Each
1330969	Large 1L Capacity 50 x 60cm Dark Green	Each




### Drycare Laminated Linen Protector

- Features non-slip waterproof backing

1290916	1L Waterproof Soft Pink 85 x 90cm	Each
---------	-----------------------------------	------

### Plastic Sheeting



1094271	Hospital sheeting PVC white 50x1m	Roll
1253752	Plastic waterproof sheet 90x120cm	Pack 100

## LAUNDRY BAGS

### Reusable Cloth Laundry Bags

- Long-lasting heavy duty polyester twill



Standard			
1274752	Blue 75 x 70cm		Each
1274805	Green 75 x 70cm		Each
1276407	Purple 75 x 70cm		Each
1274761	Red 75 x 70cm		Each
1274770	White 75 x 70cm		Each
1274799	Yellow 75 x 70cm		Each
Cytotoxic			
1281208	Cytotoxic Logo Purple 75 x 70cm		Each

### Singlet Bags

1420410	Large 520 x 280mm		Pack 100
---------	-------------------	--	----------

### Short Laundry Bags

1335561	White 75x50		Each
1348965	Blue 75x50		Each



## LINEN SKIPS

### Mobile Linen Trolleys

- Compact stainless steel tubular frame
- Foot operated lid for ease of use
- Durable acrylic lid

1323237	Single Skip Trolley
1323253	Double Skip Trolley
1323261	Triple Skip Trolley
1434373	Clean Linen Trolley Front Cover
1434403	Clean Linen Trolley Large 3 Shelf



1323253



1323253



1434373 & 1434403

### Colour ID Strips for Linen Trolleys

- Allows you to colour code linen trolleys for infection control

1371118	Blue		Each
1371126	Green		Each
1411551	Red		Each
1411578	Yellow		Each



### Impermeable Laundry Bags

- PVC coated polyester
- Waterproof and leakproof for greater infection control
- No need for pre-bagging – bag and contents are washed together on foul cycle
- Effective for approximately 200 washes
- Drying temperature is 75°C



1274500	Red 75 x 85cm		Each
1274529	Yellow 75 x 85cm		Each

### Limiting Sling

- Designed to prevent overfilling
- Results in better washing practices
- Will last for over 500 washes



1278215	Blue		Each
1277146	Green		Each
1295039	Purple		Each
1277138	Red		Each
1276840	White		Each
1278207	Yellow		Each

### Soluble Bag Liners

- Designed to fit into most laundry bags as an impermeable liner
- Manufactured from polyethylene

Soluble Seam Liners			
1274485	Red 72 x 99cm		Ctn 250
1274494	Yellow 72 x 99cm		Ctn 250
1408054	Purple 72 x 99cm		Ctn 250
Fully Soluble Liners			
1286501	Clear Aquafilm 660 x 838mm		Ctn 200




**Advise colours required**

1374710 Custom Spray Paint for Linen Trolley Lid



## SYRINGE PUMPS



### BD BodyGuard T Syringe Pump

The BD BodyGuard™ T syringe pump is a small, lightweight ambulatory infusion pump which offers flexibility and reliability. Effortless transition between hospital and home care requirements is the core of the BD BodyGuard™ T syringe pump.

1453157	BD BodyGuard T Syringe Pump
1451316	9V Battery Model 6LR61 for BodyGuard T
1264077	Lockbox for BodyGuard T
1330381	Label Subcutaneous Use 60mm x 50mm Brown
1309021	Label Line Subcutaneous Brown
1015349	I.V. Admin Extension Tubing 150cm Luer Lock



### Nipro Surefuser 100 ML X 2 DAY 2.1 ML / HR Box 10

Simple, ergonomic design makes the Surefuser+ easy to fill, easy to connect, and even easier for your patient to manage. Simply fill Surefuser+ with a standard 50ml Luer Lock syringe, attach Surefuser to IV access and open the clamp.

- Wide variety of infusion times from 30mins to 7 days
- Variety of flow rates
- Connects easily to an IV line via standard Luer connector
- Integrated line with flow regulator and 0.2µ self-venting filter
- Colour-coded for easy identification

1430610	Nipro Surefuser 100 ML X 2 DAY 2.1 ML / HR Box 10
---------	---

## MEDICATION ROUNDS

### Silent Knight Mechanical Pill Crusher

- Simple to operate – reduces repetitive injury
- Pouches contain tablet during crushing – no mess
- Portable and sturdy – easy to store on medication trolleys

1240914	Silent Knight Precision Pill Crushing Machine	Each
1240941	Silent Knight Pill Crusher Pouches Deluxe (PC1000)	Pk 1000



### Pill Cutter & Counter

- Built in magnifier that helps you align your pill or tablet
- Safety shield completely covers the cutting blade when medication is split

1365088	Pill Cutter Ezy Cut with Safety Shield & Magnifier
1368052	Triangular Pill Counter





## Utility Pouch

1435337	Nurses Utility Pouch Medium 16x20cm Violet
1435345	Nurses Utility Pouch Medium 16x20cm Black
1435353	Nurses Utility Pouch Medium 16x20cm Pink
1435418	Nurses Utility Pouch Medium 16x20cm Blue



## Porcelain Mortar & Pestle

1072336	Mortar and Pestle 100ml 100mm
---------	-------------------------------



## Do Not Disturb Apron

1303996	Medication Apron
---------	------------------

# DIAGNOSTIC MEASURE

## BLOOD PRESSURE - AUTOMATIC



### Omron HEM 7156 with Med/Lge Cuff

- Intellisense Technology:
- Automatically applies the right amount of pressure for fast, accurate and more comfortable measurements
- 60 Data Memory with date and time
- 360 degree Intelliwrap cuff (Medium to Large Size cuff 22-42cm) Snaps easily onto your arm allowing accurate readings
- Cuff wrapping guide indicator confirms the cuff has been wrapped correctly on the users arm.

#### 5 year warranty

1433377	Omron Blood Pressure Monitor & Med/Lge Cuff HEM7156
---------	---



### Omron HBP-1320 Professional Monitor incl 2 Cuff/Bldr

- With a rugged design and a practical carrying handle , HBP-1320 is designed for use in professional setting and easily carried in everyday use. The device has been clinically proven to deliver fast and reliable results.
- Designed for routine blood pressure measurements. Enhanced usability for nurses carrying the device from bed to bed and continuous blood pressure measurements. STURDY
- Sturdy enough to survive accidental falls, and allow you to carry with ease.

#### 5 year warranty

1424335	Omron HBP-1320 Professional BP Monitor incl 2 Cuff/Bldr
---------	---



### Welch Allyn ProBP 2000 Digital BP Device

- Blood pressure readings are captured on inflation in about 20 seconds
- Meets AAMI standards for blood pressure accuracy
- Compatible with Welch Allyn FlexiPort Reusable Blood Pressure Cuffs, from child size 9 up to thigh size 13 (15 - 55 cm range)
- Device includes Adult size 11 FlexiPort Reusable Blood Pressure cuff

1407821	W.A ProBP 2000 Digital BP Device & Flexiport Cuff Adult 11
---------	--

## BLOOD PRESSURE - MANUAL



### Heine Gamma G5 Sphyg

- Latex-free.
- Exclusive pressure-distribution system.
- Pressure-proof manometer. Ensures durability and reliability.
- Micro filter. Protects valve and manometer.

1224901	Gamma G5
---------	----------



### Heine Gamma G7 Sphyg

- Shockproof. Latex-Free. Upper housing made of aluminium alloy with a shock absorbing ring protecting indicator and glass.
- Complete with adult cuff, in zipper pouch
- Large Insufflation Bulb. 50 % larger bulb volume for fast inflation of the cuff.

1224910	Gamma G7
---------	----------



### Welch Allyn Durashock DS66

- Gear-free DuraShock™ technology
- Certified accuracy to +/- 3 mmHg
- Unsurpassed reliability with a unique trigger-release design
- Withstands a 30" drop

**10 year Warranty**

1146251	Durashock DS66
---------	----------------



### Spirit Aneroid Sphyg

- Large Dial for easy reading of measurements
- Latex free
- No pin stop calibrated scale

1178845	Latex Free Spirit
---------	-------------------



### Welch Allyn 767

- Lifetime calibration warranty reduces cost of ownership

**10 year Warranty**

1222700	Sphyg 767 Stand
---------	-----------------



### Sphyg Complete Aneroid

- Measurement Range: 0-300mmHg
- Cuff Material: Nylon \* Scale Glow-in-Dark \*

1036545	SM-200
---------	--------



### Welch Allyn Flexiport Cuff Range

- Affordable and easy to use, Welch Allyn FlexiPort Reusable Blood Pressure Cuffs offer all the advantages of FlexiPort technology—including simplified connectivity—in a dependable, durable design.

1286110	Size 11 25-34cm
1287281	Size 10 20-26cm
1275689	Size 11 Long 25-34cm

1286129	Cuff Size 12 32-43cm
1274171	Long Cuff Size 12 32-43cm
1274162	Thigh Size 13 40-55cm

### Welch Allyn EcoCuff

Traditional single-patient-use cuff solutions can be expensive, wasteful, and difficult to implement. EcoCuff design gives you the benefits of switching to single-patient-use, including reduced cross-contamination while controlling costs and the volume of hospital waste.



1438468	50 each of Child, Small Adult, Adult and Large Adult Cuffs
---------	--

1439367	Size 10 Small Adult 20-28 cm
1439359	Size 11 Adult 27-38cm

## VITAL SIGNS MONITOR

### Welch Allyn Connex Spot Vital Signs Monitor

- Easy-to-use, configurable 7" colour touchscreen display
- Provides accurate vital signs measurement including blood pressure averaging, spot checking and interval monitoring
- Connects wirelessly to your EMR and sends vitals to the patient chart from the point of care to help improve efficiency
- Includes W.A Flexiport Cuff Adult Size 11

1347101	Connex Spot Vital Signs Monitor NIBP/Nonin SpO2/Pro6000
1348361	Connex Spot Vital Signs Mobile Basket Stand 7000-MS3MS3
1396501	Connex Spot Vital Signs Mobile Work Stand 7000-MWS

#### Measures:

- Blood Pressure - SureBP®,
- Pulse Oximetry - Nonin®,
- Tympanic Thermometry - Thermoscan® PRO 6000

1287281	Flexiport Cuff Small Adult Size 10 20-26cm
1286110	Flexiport Cuff Adult Size 11 25-34cm
1274171	Flexiport Lge Adult Long Cuff Size 12 32-43cm No Tub
1286129	Flexiport Large Adult Cuff Size 12 32-43cm No Tube
1274162	Flexiport Adult Thigh Size 13 40-55cm
1275689	Flexiport Cuff Adult Size 11 Long 25-34cm



## POINT OF CARE TESTING Clinitek Status +

Simple, point-of-care urinalysis with enhanced clinical information

- Easy-to-use, automated analyzer designed to deliver quick, accurate, reliable results
- Comprehensive urine test menu
- Printed record
- Checks integrity of each test strip
- Identifies type of strip automatically



1296175	Clinitek Status & Urine Chemistry Analyser
1163181	Thermal Paper Rolls for Clinitek Analyser
1312502	DCA/ Clinitek Printer Labels
1005899	MultistiX-10SG Strips 100



## PULSE OXIMETER



### Finger Pulse Oximeter

1319124	Aeon Pulse Oximeter Finger Adult LCD Multi Display
---------	--



### Heart Sure Fingertip Pulse Oximeter

1339931	Heart Sure Fingertip Pulse Oximeter
---------	-------------------------------------



## DIAGNOSTIC SETS

### Heine mini 3000 Diagnostic Set



Ideally suited to the medical professional who values maximum quality and minimum size. Conveniently compact, lightweight and with a unique on/off switch, the HEINE mini3000 series makes an ideal portable set. The pocket-sized set, still in HEINE quality

1263688

Heine mini 3000 diagnostic set

### Welch Allyn Macroview & Coaxial Set



An innovative diagnostic set that features a 3.5V Coaxial Ophthalmoscope, 3.5V MacroView Otoscope and Lithium Ion Power Handle and LED globes.

Including Desk Pod and Charger Hard carry case included 1 year warranty\*

1364960

Welch Allyn Macroview & Coaxial Set

### Welch Allyn Pocket LED Diagnostic Sets

The Welch Allyn Pocket LED diagnostic set combines function and style whilst providing maximum performance and durability in a lightweight, compact set. It is the first pocket size LED diagnostic set that uses advanced SureColour LED technology, which delivers true tissue colour rendition.



1342738

Set 2 Handles in Soft Case Blue

1342703

Set 2 Handles in Soft Case Black

1412159

Otoscope Only Blue No Case

1342967

Otoscope Only in Soft Case Blue

## THERMOMETERS

### Welch Allyn Braun Thermoscan Pro 6000

- Accurate, reliable and comfortable temperature readings
- Innovative PerfectTemp technology adjusts for variability in probe placement
- ExacTemp technology detects the stability of the probe during measurement
- 60 second pulse timer assists with the manual measurement of pulse rate, respiration

1337564	W.A Thermoscan PRO6000 with Small Cradle
1337580	W.A Thermoscan PRO6000 Charging Dock
1337572	W.A Rechargeable Battery Pack for Thermoscan PRO6000
1337599	W.A Thermoscan Probe Covers to Suit PRO6000/PRO4000



### Omron TH839S Ear Thermometer

- 1 second quick measurement
- Fever alarm with beeper for body temperature
- 9 test memory
- Ergonomic design for easy handling

1307495	Omron Ear Thermometer Gentle Temp TH839S
1307754	Thermometer Ear Probe Covers Omron TH839S



### Infrared Forehead Thermometer

- Measure Distance is between 1-3cm
- Hold 20 readings in memory.
- Colour coded temperature reading.

1431803	Infrared Forehead Thermometer
---------	-------------------------------



## BLADDER SCANNERS



### Portable Cubescan BioCon-900

- The BioCon-900 is the latest release from Mcube, offering a wireless, non-invasive 3D ultrasound system that measures the volume of urine in a patient's bladder.
- The BioCon-900 uses ultrasound signals to determine the bladder's outline and total volume of urine. The revolutionary real-time Pre-Scan function locates the bladder prior to scanning and improves the accuracy of results.
- Immediate printout of results via the internal printer or use the included CubeScan PC software to upload saved data to a computer for later review. 12 Month Warranty

1406426	BioCon 900D Bladder Scanner with docking station
---------	--



### Bladder Scanner Touch Screen Printer Kaixin BVT01

- The portable bladder scanner BVT01 is a non-invasive ultrasound scanner used to measure bladder volume for assessing urinary retention and urinary incontinence. It aims to reduce the conduction time urinary frequency and reduce the risk of urinary tract infection and for treatment of urinary system diseases.
- Features smart and compact design with touchscreen, simple operation and advanced calculating method ensures the correct result for every anomaly shape.
- 24 Month Warranty main unit
- 12 Months Warranty Probes

1428101	Bladder Scanner + Touch Screen + Printer Kaixin BVT01
---------	---

## STETHOSCOPES

### Spirit Dual Head

- Grandeur series lightweight nurse single head stethoscope
- Features chest piece made of aluminium material with metal retaining ring, simplex; binaural and 1 piece molded PVC tubing

1203198	Stethoscope Nurses Complete Orange Tube Spirit
1203385	Stethoscope Spirit Deluxe Dual Head Neon Orange
1153069	Stethoscope Spirit Deluxe Dual Head Pink
1153078	Stethoscope Spirit Deluxe Dual Head Aqua Blue
1153031	Stethoscope Spirit Deluxe Dual Head Red
1153040	Stethoscope Spirit Deluxe Dual Head Burgundy



### Liberty Nurses Single Head

- Grandeur series lightweight nurse single head stethoscope
- Features chest piece made of aluminium material with metal retaining ring, simplex; binaural and 1 piece molded PVC tubing

1314009	Dual Head Complete Black Liberty	1043049	Single Head Royal Blue Liberty
1316540	Dual Head Navy Liberty	1043058	Single Head Teal Liberty
1316486	Dual Head Teal Liberty	1043076	Single Head Red Liberty
1316508	Dual Head Purple Liberty	1313932	Single Head Black Liberty
1316494	Dual Head Magenta Liberty	1324675	Single Head Purple Liberty
1316540	Dual Head Navy Liberty		
1316508	Dual Head Purple Liberty		



## SCALES



### Healthweigh Wheelchair Scale

The Healthweigh® Wheelchair Scale allows practices to measure patients confined to a wheelchair in a safe and comfortable manner. Patients can be placed onto the scale via the scale's two wide accessibility ramp. The Healthweigh® Wheelchair Scale also includes movement compensation technology, which eliminates involuntary movement during the weighing process to ensure accurate results. The scale also features a manual/automatic tare option, which allows the weight of wheelchair to be reduced in calculations.

Code	Capacity	Size	Graduation
1392719	360kg	1250x940x82mm	0.1kg



### Weighing Chair

- Dry battery or AC power
- Comfortable white plastic seat
- Full digital calibration (FDC) ensuring accurate and simple calibration
- Auto power off to preserve battery life
- Safety wheels with brakes, standard foot rest

Code	Capacity
1049786	150KG

### Weighing Chair Scale

Code	Capacity
1360914	300kg



### Seca Weighing Chair 952

Code	Capacity
1330691	200KG



## EXAM LIGHTS

### Welch Allyn GS 300 General Exam Light

- Single White LED
- 20,000 Lux at typical working distance
- Perfect combination of intense light output and cool operation
- Touchless on/off
- Ergonomic power connection
- Affordable price and rugged performance

1280554	W.A LED General Examination Light with Mobile Stand GS300
1280562	W.A LED General Examination Light with Table/Wall Mnt GS300
1297732	W.A Extendable Wall Mount for GS Lights



## STAND ALONE ECG

### Welch Allyn® CP 150

Stand Alone ECG with Interpretation

- 7" touch-screen display and keyboard
- Advanced filters for optimal ECG trace quality
- Three lead ECG preview
- Battery or mains operated for ultimate flexibility

1314270	W.A CP150 Interpretive 12 Lead Resting ECG Touch Screen
1316516	W.A CP150 Hospital Cart W/Out Cable Arm/Shelf
1316524	W.A Cable Arm & Shelf Option for CP150 Cart



### ANAEMIA SCREENING HemoCue HB301 Analyser

- Provides measured haemoglobin within 3 seconds
- Small sample size - 10µl
- Acceptable samples: capillary, venous or arterial
- Consumables stable at room temperature

1397931	Hemocue HB301 Haemoglobin Analyser
1379984	Hemocue HB301 Microcuvettes (4 x 50)



### GLUCOMETER Accu-Chek Instant S

- Intuitive 9 colour target range indicator helps understand results easily
- Proven accuracy
- 1 button navigation
- Easy-edge blood dosing, no set up

1410105	Accu-Chek Instant S Glucometer Meter
1400622	Accu-Chek Instant S Blood Glucose Strips B50
1410202	Accu-Chek Instant Controls Solution INTL II 2 x 2.5ml



### COMPACT 525 - 5 LITRE OXYGEN CONCENTRATOR

Delivering up to 5 litres per minute with high oxygen concentration across all flow rates and adopting the latest developments in the reduction of sound and power consumption, the latest Compact 525 increases patient comfort and lowers operating costs.

- Recessed humidifier nook prevents damage
- New maintenance free innovative rotary valve
- Intelligent power management system utilises turn-down® technology providing less power consumption below 2 l/min flow rates
- Weighing only 16.3 kg the Compact 525 is the smallest and lightest 5 litre concentrator from DeVilbiss Healthcare
- Circuit breaker with reset function
- Oxygen Sensing Device for continuously monitoring Oxygen levels
- Visual and audible alarms for low Oxygen levels, power failure, pressure drop and service required

1442759	Compact 525 - 5 Litre Oxygen Concentrator
---------	---



## SUCTION PUMPS

### Askir 30 Surgical Aspirator/ Askir 36BR

- Easily transportable
- Highly heat resistant and electrically insulated plastic
- Overflow valve and vacuum indicator

1272963	Askir Suction Pump 30 Ltr 230V (1 Ltr Jar)
1323822	Askir 36BR Suction Pump with Built in Batter Pack
1276354	Askir Antibacterial Filter (also fits Liberty)
1309595	Askir 1000ml Jar (Autoclavable)



### Disposable Items for Askir 30/36BR

1439952	Askir Round Spacer For Suction Unit
1323814	Flovac Askir Reusable Container 1Lt (Takes Liner 1313150)
1313150	Flovac Disposable Liners fits Askir 1ltr Canister
1439944	Flovac Conical Connector (Yellow) for Flovac Disposable Lid

### Replacement Items For Askir 30/36BR

1276354	Askir Antibacterial Filter (also fits Liberty)
1428441	Askir Suction Tubing Set 3 Tubes 6x10x140cm(1) 6x10x13cm(2)



## NEBULISER THERAPY

### AiroMax Nebuliser Therapy Unit

- Nebulisation rate 0.55ml/minute
- Five year warranty
- Compressor max pressure 2.5 bar, Power 240V, 50Hz, 1.0 Amp
- Noise level 52 dbA
- Masks, bowl, tubing, filters and carry case
- Quiet piston operation
- Allows patient-specific individualisation



1152348	AiroMax Nebuliser Therapy Unit
1417177	Filter for Airo MAX Nebuliser Therapy Unit PK3

## DEFIBRILLATOR

### LIFEPAK CR2 ESSENTIAL SEMI-AUTOMATIC

- Semi-automatic configuration
- Compact and easy to use lifesaving technology
- Easy to follow clear, calm verbal guidance
- Easy two step operation
- Includes two sets of adult electrode pads & first responder kit
- Eight year manufacturers warranty



1415743	Defib Lifepak CR2 Essential Semi-Automatic
1419005	CR2 Essential Lithium Replacement Battery
1419013	CR2 Essential Wall Mount Bracket - White
1418998	CR2 Quik-Step Essential Replacement Electrodes



**SSS AUSTRALIA**  
HEALTHCARE SUPPLIES

 **1800 777 518**

 **sales@sssaustralia.com.au**

 **www.sssaustralia.com.au**

All orders may be subject to a shipping fee. Pictures for illustration purposes only. All items while stocks last. E&O.E. Additional 10% GST will be added to all prices shown, unless marked with an \* (asterisk denotes GST Free). Pricing listed RRP and indicative only.