



AGANI™ Hypodermic Needles

Smooth and easy puncture for patient comfort.

AGANI™ Hypodermic Needles

LATEX FREE

Easy to use. Individually packaged in a sterile colour-coded blister pack for easy gauge identification.

Suitable for all applications. Wide range of sizes (18G to 30G).

Additional features. Visual observation of flashback through semi-transparent colour-coded needle hub. Thin wall needles allow you to choose a smaller gauge size with increased flow rates over regular walled needles.

Agani™ Hypodermic Needles

Product Code	Colour	Gauge	Length	Carton
			mm	
AN*1838R1		18G	38	10000
AN*1925R1		19G	25	10000
AN*1938R1		19G	38	10000
AN*2025R1		20G	25	10000
AN*2038R1		20G	38	10000
AN*2116R1		21G	16	10000
AN*2125R1		21G	25	10000
AN*2138R1		21G	38	10000
AN*2219R1		22G	19	10000
AN*2225R1		22G	25	10000
AN*2232R1		22G	32	10000
AN*2238R1		22G	38	10000
AN*2316R1		23G	16	10000
AN*2325R1		23G	25	10000
AN*2332R1		23G	32	10000
AN*2425R1		24G	25	10000
AN*2516R1		25G	16	10000
AN*2525R1		25G	25	10000
AN*2613R1		26G	13	10000
AN*2623R1		26G	23	10000
AN*2716R1		27G	16	10000
AN*2719R1		27G	19	10000
AN*3013R1		30G	13	10000

Shelf Pack = 100

Needle Length: 50mm

Product Code	Colour	Gauge	Length	Carton
			mm	
AN*1850R1		18G	50	5000
AN*1950R1		19G	50	5000
AN*2050R1		20G	50	5000
AN*2150R1		21G	50	5000
AN*2250R1		22G	50	5000

Shelf Pack = 50

Drawing Up Needle (Blunt Needle)

AN*1838N		18G	38	10000
----------	--	-----	----	-------

Shelf Pack = 100

Terumo Australia Pty Ltd

Macquarie Park NSW 2113 Australia
 Mount Albert Auckland 1025 New Zealand
 E: cs_australia@terumo.co.jp

T: 1800 837 866
 T: 0800 66 77 57

AGANI™ trademark of Terumo Corporation.
 © 2018 Terumo Australia Pty Ltd. GHP180223D651V5 | PBD00276

