



## SSS Australia Fitment Guide

The SureGard® filter is an Australian owned, highly effective filter barrier device. The SureGard® provides both technicians and patients with optimal protection from cross contamination and infection. The R&D of the SureGard® is attributed to contributions from the medical and scientific community of Australia and New Zealand.

### Purple RJVKB2 B25 #1437151

Ex. Diameter: 33mm-32mm  
Int. Diameter: 29mm

- Fits MIR Disposable (with cardboard mouthpiece piece attached)  
Code: #1342673
- Fits QRS Orbit Turbine  
Code : #1262804

### Green RJVKB4 B20 #1204524 B25 #1434365

Ex. Diameter: 32.5mm  
Int. Diameter: 29mm

- Fits MIR Reuseable Turbine  
Code: #1217658

### Pink RJVKB1 B20 #1205192 B25 #1434357

Ex. Diameter: 28.5mm  
Int. Diameter: 25mm

- Fits QRS Orbit (tight fit)  
Code: #1262804
- Fits MIR Reuseable Turbine  
Code: #1217658

### Blue RJVKB3 B25 #1435159

Ex. Diameter: R:33mm-32mm L:28.5mm  
Int. Diameter: R:29mm L:25mm

- Fits MIR Reuseable Turbine  
Code: #1217658
- Fits QRS Orbit  
Code: #1262804

### Pearl RJVKB7 By Request Only

Ex. Diameter: 20mm  
Int. Diameter: 18.5mm

- Fits EasyOne Spirometer (spirettes)  
- not currently stocked by SSS Australia.

### Red RJVKB5 By Request Only

Ex. Diameter: 31mm  
Int. Diameter: 29mm

- Fits MediSoft Systems - not currently stocked by SSS Australia

### White RJVKB9 By Request Only

Ex. Diameter: 30mm  
Int. Diameter: 28mm

- Fits MIR Reuseable Turbine (Tight Fit)
- Fits MIR Disposable with Cardboard removed
- CosMed Pony, CosMed MicroQuark - not currently stocked by SSS Australia

### Orange RJVKB8 By Request Only

Ex. Diameter: 36mm  
Int. Diameter: 34mm

- MedGraphics systems - not currently stocked by SSS Australia.

### Clear RJVKB6 B100 #1230177

Ex. Diameter: R:25.5mm L:22mm  
Int. Diameter: R:22mm L:19.5mm

- Ventilation systems, nebulisers, CPAP & BiPAP equipment.

**P: (07) 3890 2600**

**www.sssaustalia.com.au**  
**sales@sssaustalia.com.au**