

∴ HEINE® NC 1 Dermatoscope

Non-contact examination with unrivalled quality.



OPTIONAL
CONTACT PLATE



POLARISED
LIGHT



LED NOW IN HEINE QUALITY.



HEINE® NC 1 DERMATOSCOPE. THE BEST OF BOTH WORLDS.

The new NC1 Dermatoscope is truly two instruments in one. It allows detailed, non-contact examination of the patient in a first screening, without unnecessary skin contact and disinfection. Should a skin lesion need to be examined more closely, a contact plate can be added, transforming the NC 1 into a contact Dermatoscope with increased magnification and a scale. The high quality aspherical optics provides unmatched sharpness and resolution; LED^{HQ} illumination offers true and natural colour rendering. The solid aluminum frame combined with HEINE Quality guarantee a long lifespan.



NON-CONTACT DERMATOSCOPY

Non-contact screening for efficient examination and a professional distance from the patient. Large field of view, excellent depth of field with 6x magnification.



POLARISED LIGHT

Detailed views of crystalline and vascular structures without reflexes.



OPTIONAL CONTACT PLATE

Magnetic contact plate with 9x magnification provides high-quality contact images for a more detailed diagnosis.

LED IN HEINE QUALITY

Advanced LED illumination for an absolutely homogeneous light with true and natural colour rendering.



NC1 Dermatoscope

NC1 Dermatoscope Head, without handle	D-008.78.130
NC1 contact plate with scale	D-000.78.115

NC1 Dermatoscope Sets

NC1 Dermatoscope Head on mini3000 battery handle with batteries	D-008.78.132
NC1 Dermatoscope Head on mini3000 rechargeable handle with rechargeable battery	D-008.78.133
NC1 Dermatoscope Head on mini3000 battery handle with contact plate and batteries	D-008.78.134
NC1 Dermatoscope Head on mini3000 rechargeable handle with contact plate and rechargeable battery	D-008.78.135
NC1 Dermatoscope Head, contact plate with scale, mini3000 LED F.O. Otoscope Head, 2 mini3000 rechargeable handles with rechargeable battery, mini NT Table Charger [01]	D-890.11.022
NC1 Dermatoscope Head without contact plate, mini3000 LED F.O. Otoscope Head, 2 mini3000 battery handles with battery and 5 of both 2.5 and 5mm dia. AllSpec disposable tips, Zipper case [02]	D-891.11.021



without contact plate



with contact plate



[01]



[02]

