

The Haemolance Plus®

Product range

Haemolance Plus® Micro Flow

The Micro Flow is ideal for diagnostic testing requiring a sample less than four microliters of blood. The penetration depth is 1.6 mm with a 0, 4 mm 28-gauge needle. Typical test include blood glucose.

Product #: 53159 - Box of 150, 53158 - Box of 50



Haemolance Plus® Low Flow

The Low Flow is ideal for diagnostic testing requiring a small sample size of blood but more than four microliters. The penetration depth is 1.4 mm with a 25-gauge needle.

Product #: 53156 - Box of 150, 53155 - Box of 50



Haemolance Plus® Normal Flow

The Normal Flow is ideal for diagnostic testing requiring a larger blood sample. The penetration depth is 1.8 mm with a 21-gauge needle. Typical tests include hemoglobin and cholesterol testing.

Product #: 53152 - Box of 150, 53151 - Box of 50



Haemolance Plus® High Flow

The High Flow is ideal for blood banks, labs or places that want greater than 10 microliters of capillary blood. The penetration depth is 1.8 mm with a 18-gauge needle.

Product #: 51361 - Box of 150, 51362 - Box of 50



Haemolance Plus® Max Flow

This blade style lancet is ideal for blood banks and laboratories or places where you need a large amount of blood. The blade has a penetration depth of 1.6 mm.

Product #: 53169 - Box of 150



Haemolance Plus® Pediatric

This blade style lancet is ideal for pediatric heel sticks on neonates. The blade has a penetration depth of 1.2 mm.

Product #: 53164 - Box of 150

