

GAMMEX[®] Latex Powdered

S U R G I C A L G L O V E S

Latex, powdered surgical glove delivering comfort and barrier protection

KEY FEATURES AND BENEFITS:

- Uses Ansell's exclusive glove formulation delivering strong, reliable gloves
- Textured finish ensures a secure grip for handling instruments under dry and damp conditions
- Made with Ansell's unique anatomic formers to enhance comfort and fit to help reduce hand fatigue
- Powdered for ease of donning

Ansell recommends the use of powder-free gloves to minimize adverse effects for both healthcare workers and patients.

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.



GAMMEX[®] Innovation that Protects

Driven by innovation, GAMMEX[®] surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX[®] Latex Powdered

S U R G I C A L G L O V E S

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powdered
Color	Natural
Thickness	Standard
Cuff Style	Straight
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Powdered
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	No, not tested in accordance with ASTM D6978
Freedom from Holes (Inspection level I)	0.65 AQL
Product Standards	AS/NZS 4179, EN 420:2003 + A1:2009, EN 455 1-4, ISO 10282, JIS T9107
Quality / Environmental Standards	EN556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	<ul style="list-style-type: none"> Zinc Diethyldithiocarbamate (ZDEC) Zinc 2-mercaptobenzothiazole (ZMBT)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Ansell Product Code	330047055	330047060	330047065	330047070	330047075	330047080	330047085	330047090

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 1-844-494-7854
 Email: insidesalesus@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu

Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21)2407 2286/88/62
 Email: infochina@ansell.com