



## Need a Sling?

Find the correct sling for your manual handling scenario on page 16



## Quick Tech Specs



See page 14 for detailed specs

## PATIENT LIFTER - STEEL

A robust workhorse, the Aidacare Aspire 200 is designed to meet the rigorous demands of professional and facility environments.



### ERGONOMIC HANDLES

Provide multiple grip points to suit carers of different heights and strengths



### ERGONOMIC HANDSET

The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling



### POWERED LEG FUNCTION

Opens and closes legs for lifting stability or accessibility



### QUICK RELEASE YOKE ATTACHMENT

One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety



### BUILD QUALITY

Strong powder coated steel frame with industry leading Linak electronics and warranty provision



### SAFETY

Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart



### INTELLIGENT DIAGNOSTICS

State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance



### QUALITY CASTORS

2 x 75mm front castors for maximum under bed clearance and accessibility.  
2 x 100mm rear locking castors for safety and manoeuvrability



PRODUCT CODE: L5S390620